

CARETAKER LIBRARY API 1.5.3 REFERENCE MANUAL

Caretaker Medical

May 25, 2018

Contents

1	Introduction	1
2	Library Integration	3
3	Monitoring the Device	7
4	Module Index	13
4.1	Modules	13
5	Data Structure Index	15
5.1	Data Structures	15
6	Module Documentation	17
6.1	Device Information	17
6.1.1	Detailed Description	18
6.1.2	Macro Definition Documentation	18
6.1.2.1	libct_device_get_state	18
6.1.2.2	libct_device_uninitialized	18
6.1.2.3	libct_device_discovering	19
6.1.2.4	libct_device_connecting	19
6.1.2.5	libct_device_connected	19
6.1.2.6	libct_device_disconnecting	19
6.1.2.7	libct_device_disconnected	20
6.1.2.8	libct_device_monitoring	20
6.1.2.9	libct_device_measuring	20
6.1.2.10	libct_device_get_class	21

6.1.2.11	libct_device_get_name	21
6.1.2.12	libct_device_get_address	21
6.1.2.13	libct_device_get_serial_number	21
6.1.2.14	libct_device_get_hw_version	22
6.1.2.15	libct_device_get_fw_version	22
6.1.2.16	libct_device_get_context	22
6.2	Primary API	23
6.2.1	Detailed Description	23
6.2.2	Function Documentation	23
6.2.2.1	libct_init()	23
6.2.2.2	libct_deinit()	24
6.2.2.3	libct_start_discovery()	24
6.2.2.4	libct_stop_discovery()	25
6.2.2.5	libct_connect()	25
6.2.2.6	libct_disconnect()	26
6.2.2.7	libct_start_monitoring()	27
6.2.2.8	libct_stop_monitoring()	27
6.2.2.9	libct_start_measuring()	28
6.2.2.10	libct_stop_measuring()	28
6.3	Secondary API	29
6.3.1	Detailed Description	30
6.3.2	Macro Definition Documentation	30
6.3.2.1	libct_dp_count	30
6.3.2.2	libct_get_dp	31
6.3.2.3	libct_get_last_dp	31
6.3.2.4	libct_get_first_dp	32
6.3.2.5	for_each_dp	32
6.3.2.6	libct_inc_cuff_pressure	32
6.3.2.7	libct_dec_cuff_pressure	33
6.3.3	Function Documentation	33

6.3.3.1	libct_get_device()	33
6.3.3.2	libct_set_app_specific_data()	34
6.3.3.3	libct_get_app_specific_data()	34
6.3.3.4	libct_get_version_string()	35
6.3.3.5	libct_get_build_date_string()	35
6.3.3.6	libct_set_log_level()	35
6.3.3.7	libct_recalibrate()	36
6.3.3.8	libct_adjust_cuff_pressure()	36
6.3.3.9	libct_rd_cuff_pressure()	37
6.3.3.10	libct_vent_cuff()	37
6.3.3.11	libct_clr_status()	38
6.3.3.12	libct_diag_flush()	38
6.3.3.13	libct_wrt_snr_min()	38
6.3.3.14	libct_rd_snr_min()	39
6.3.3.15	libct_wrt_display_state()	39
6.3.3.16	libct_rd_display_state()	40
6.3.3.17	libct_wrt_recal_itvl()	40
6.3.3.18	libct_rd_recal_itvl()	40
6.3.3.19	libct_wrt_waveform_clamping()	41
6.3.3.20	libct_rd_waveform_clamping()	41
6.3.3.21	libct_rd_vitals_filter()	42
6.3.3.22	libct_wrt_vitals_filter()	42
6.3.3.23	libct_wrt_simulation_mode()	42
6.3.3.24	libct_wrt_motion_timeout()	43
6.3.3.25	libct_rd_motion_timeout()	43
6.3.3.26	libct_rd_persistent_log()	44

7 Data Structure Documentation	45
7.1 libct_app_callbacks_t Struct Reference	45
7.1.1 Detailed Description	46
7.1.2 Field Documentation	47
7.1.2.1 on_device_discovered	47
7.1.2.2 on_discovery_timedout	48
7.1.2.3 on_discovery_failed	48
7.1.2.4 on_device_connected_not_ready	49
7.1.2.5 on_device_connected_ready	50
7.1.2.6 on_connect_error	50
7.1.2.7 on_connect_timedout	51
7.1.2.8 on_device_disconnected	52
7.1.2.9 on_start_monitoring	52
7.1.2.10 on_stop_monitoring	53
7.1.2.11 on_start_measuring	54
7.1.2.12 on_stop_measuring	54
7.1.2.13 on_data_received	55
7.1.2.14 on_data_error	56
7.1.2.15 on_rd_snr_min_rsp	56
7.1.2.16 on_wrt_snr_min_rsp	57
7.1.2.17 on_rd_display_state_rsp	57
7.1.2.18 on_wrt_display_state_rsp	58
7.1.2.19 on_rd_recal_itvl_rsp	59
7.1.2.20 on_wrt_recal_itvl_rsp	59
7.1.2.21 on_rd_cuff_pressure_rsp	60
7.1.2.22 on_vent_cuff_rsp	61
7.1.2.23 on_clr_status_rsp	61
7.1.2.24 on_diag_flush_rsp	62
7.1.2.25 on_wrt_waveform_clamping	62
7.1.2.26 on_rd_waveform_clamping	63

7.1.2.27	on_rd_vitals_filter	64
7.1.2.28	on_wrt_vitals_filter	64
7.1.2.29	on_rd_motion_timeout	65
7.1.2.30	on_rd_persistent_log	65
7.2	libct_battery_info_t Class Reference	66
7.2.1	Detailed Description	66
7.2.2	Field Documentation	66
7.2.2.1	voltage	67
7.2.2.2	timestamp	67
7.3	libct_bp_settings_t Class Reference	67
7.3.1	Detailed Description	67
7.3.2	Field Documentation	67
7.3.2.1	systolic	67
7.3.2.2	diastolic	68
7.4	libct_cal_curve_t Class Reference	68
7.4.1	Detailed Description	68
7.4.2	Field Documentation	68
7.4.2.1	valid	69
7.4.2.2	data_id	69
7.4.2.3	val1	69
7.4.2.4	val2	69
7.4.2.5	val3	69
7.5	libct_cal_t Struct Reference	69
7.5.1	Detailed Description	70
7.5.2	Field Documentation	70
7.5.2.1	type	70
7.5.2.2	posture	70
7.5.2.3	config	70
7.6	libct_cal_type_t Class Reference	71
7.6.1	Detailed Description	71

7.7	libct_context_t Class Reference	71
7.7.1	Detailed Description	71
7.8	libct_cuff_pressure_t Class Reference	71
7.8.1	Detailed Description	72
7.8.2	Field Documentation	72
7.8.2.1	valid	72
7.8.2.2	value	72
7.8.2.3	target	72
7.8.2.4	snr	73
7.8.2.5	timestamp	73
7.9	libct_device_class_t Class Reference	73
7.9.1	Detailed Description	73
7.10	libct_device_state_t Class Reference	73
7.10.1	Detailed Description	73
7.11	libct_device_status_t Class Reference	74
7.11.1	Detailed Description	75
7.11.2	Field Documentation	75
7.11.2.1	simulation_enabled	75
7.11.2.2	inflated_indicator	76
7.11.2.3	clock_wrap_around	76
7.11.2.4	battery_voltage_low	76
7.11.2.5	pump_overrun	76
7.11.2.6	ble_temperature_sensor_paired	76
7.11.2.7	ble_stream_control	76
7.11.2.8	manual_cal_mode	77
7.11.2.9	motion_event	77
7.11.2.10	poor_signal	77
7.11.2.11	data_valid	77
7.11.2.12	calibration_offset_failed	77
7.11.2.13	weak_signal	77

7.11.2.14 bad_cuff	78
7.11.2.15 ble_adv	78
7.11.2.16 recal_soon	78
7.11.2.17 too_many_fails	78
7.12 libct_device_t Class Reference	78
7.12.1 Detailed Description	79
7.12.2 Field Documentation	79
7.12.2.1 get_state	79
7.12.2.2 get_class	80
7.12.2.3 get_name	80
7.12.2.4 get_address	81
7.12.2.5 get_serial_number	81
7.12.2.6 get_hw_version	83
7.12.2.7 get_fw_version	83
7.12.2.8 get_context	84
7.13 libct_init_data_t Struct Reference	84
7.13.1 Detailed Description	84
7.13.2 Field Documentation	85
7.13.2.1 device_class	85
7.14 libct_monitor_flags_t Class Reference	85
7.14.1 Detailed Description	85
7.15 libct_param_pulse_t Class Reference	85
7.15.1 Detailed Description	86
7.15.2 Field Documentation	86
7.15.2.1 valid	86
7.15.2.2 protocol_header	86
7.15.2.3 t0	87
7.15.2.4 t1	87
7.15.2.5 t2	87
7.15.2.6 t3	87

7.15.2.7	p0	87
7.15.2.8	p1	87
7.15.2.9	p2	88
7.15.2.10	p3	88
7.15.2.11	ibi	88
7.15.2.12	as	88
7.15.2.13	sqe	88
7.15.2.14	pressure	88
7.15.2.15	time	89
7.15.2.16	waveform_len	89
7.16	libct_posture_t Class Reference	89
7.16.1	Detailed Description	89
7.17	libct_pulse_ox_t Class Reference	89
7.17.1	Detailed Description	90
7.17.2	Field Documentation	90
7.17.2.1	sao2	90
7.17.2.2	pulse_rate	90
7.17.2.3	timestamp	90
7.18	libct_pulse_t Class Reference	90
7.18.1	Detailed Description	91
7.18.2	Field Documentation	91
7.18.2.1	timestamp	91
7.19	libct_pulse_waveform_t Struct Reference	91
7.19.1	Detailed Description	92
7.19.2	Field Documentation	92
7.19.2.1	receive_time	92
7.20	libct_status_t Class Reference	92
7.20.1	Detailed Description	92
7.21	libct_stream_data_t Class Reference	92
7.21.1	Detailed Description	94

7.21.2	Field Documentation	95
7.21.2.1	device	95
7.21.2.2	device_status	95
7.21.2.3	battery_info	95
7.21.2.4	datapoints [1/8]	95
7.21.2.5	count	95
7.21.2.6	vitals	96
7.21.2.7	datapoints [2/8]	96
7.21.2.8	cuff_pressure	96
7.21.2.9	datapoints [3/8]	97
7.21.2.10	temperature	97
7.21.2.11	datapoints [4/8]	97
7.21.2.12	pulse_ox	97
7.21.2.13	datapoints [5/8]	98
7.21.2.14	vitals2	98
7.21.2.15	datapoints [6/8]	98
7.21.2.16	raw_pulse	99
7.21.2.17	int_pulse	99
7.21.2.18	datapoints [7/8]	99
7.21.2.19	param_pulse	100
7.21.2.20	datapoints [8/8]	100
7.21.2.21	cal_curve	100
7.21.2.22	receive_time	101
7.22	libct_temperature_t Class Reference	101
7.22.1	Detailed Description	101
7.22.2	Field Documentation	101
7.22.2.1	value	101
7.22.2.2	timestamp	102
7.23	libct_version_t Struct Reference	102
7.23.1	Detailed Description	102

7.24 libct_vitals2_t Class Reference	102
7.24.1 Detailed Description	103
7.24.2 Field Documentation	103
7.24.2.1 blood_volume	103
7.24.2.2 timestamp	103
7.25 libct_vitals_t Class Reference	104
7.25.1 Detailed Description	104
7.25.2 Field Documentation	105
7.25.2.1 valid	105
7.25.2.2 bp_status	105
7.25.2.3 systolic	105
7.25.2.4 diastolic	105
7.25.2.5 map	105
7.25.2.6 heart_rate	106
7.25.2.7 respiration	106
7.25.2.8 as	106
7.25.2.9 sqe	106
7.25.2.10 timestamp	106

Chapter 1

Introduction

The Caretaker library is a cross-platform library to link with Android, Linux, and Windows applications. The library provides an interface to the Caretaker 4 wireless vitals signs monitoring device.

Scope

This version of the manual covers the Caretaker Library Linux API. The Caretaker library Android and Windows APIs are covered in separate documents. Please contact customer service for copies of the other APIs documentation as needed.

Getting Started

Read the following sections to integrate the library with an application to monitor the Caretaker device.

- [Library Integration](#)
- [Monitoring the Device](#)

Chapter 2

Library Integration

This section provides an overview of the library integration with Linux applications.

Library Package

Linux versions of the library are released in zip format *ctlibrary-linux-ARCH-VERSION.zip* for use with Linux applications, where *ARCH* is the target platform architecture, and *VERSION* is the library version bundling the following components.

Component	Description
caretaker.h	Header file exposing the Caretaker APIs to the application.
libcaretaker.a	Caretaker library image to use for static linking with the application.
libcaretaker.so	Caretaker library image to use for dynamic linking with the application.
libbluetooth.a	Precompiled version of the BlueZ Linux Bluetooth library for static linking provided for convenience.
libbluetooth.so	Precompiled version of the BLueZ Linux Bluetooth library for dynamic linking provided for convenience.
manual	The library API documentation. Documentation is available in pdf and html formats. The pdf version is supplied for completeness though the formatting is optimized for viewing in html format. For viewing the html version main page open html/index.html .
examples	Examples illustrating library use.

Supported Platforms

This library version links with Ubuntu and Debian applications running on x86 and ARM based platforms.

The library was tested on the following Linux systems.

- Ubuntu 14.04 Desktop
- Rasberry Pi running Raspbian GNU/Linux 9

Additional Hardware

If the Linux computer that will run the application does not have Bluetooth support or the Bluetooth hardware is not current, you will need the CSR8510 A10 BLE dongle or similar to enable the computer to communicate with the Caretaker device. If so, plug the dongle into a USB port on the computer and then run the following commands to configure it. Replace *hcin* with the actual host control interface name printed in the *hcitool* output for the dongle.

```
hcitool dev
sudo hciconfig <hcin> reset
```

Linking with the Library

Linux applications can link the with the library using either static or dynamic linking.

Static Linking

For static linking, the static library image *libcaretaker.a* and the supplied *libbluetooth.a* for BLE communication will need to be linked to the application. Copy these two library images to the application project library directory and link them as follows if compiling with *gcc*.

```
$(PROG): $(OBJS)
gcc -static -o $@ $^ $(LIBS)/libcaretaker-VERSION.a $(LIBS)/libbluetooth.a -lpthread
```

Where *PROG* is the application program name, *OBJS* are the application object files, and *LIBS* is the path to the project library directory. Note the caretaker library uses POSIX threads so *-lpthread* is required to link the pthread library.

Dynamic Linking

For dynamic linking, the dynamic library image *libcaretaker.so* and the supplied *libbluetooth.so* for BLE communication will need to be linked to the application. Copy these library images to the application project library directory and link them as follows if compiling with *gcc*.

```
$(PROG): $(OBJS)
gcc -o $@ $^ -lcaretaker-linux-VERSION -lbluetooth -lpthread -L$(LIBS)
```

Where *PROG* is the application program name, *OBJS* are the application object files, and *LIBS* is the path to the project library directory. Note the caretaker library uses POSIX threads so *-lpthread* is required to link the pthread library.

Debugging

The library supports generating log messages with the following levels of verbosity.

Log Level	Description
0	Show all log messages. This is the most verbose level.
1	Show informational, warning, and error messages only.
2	Show warning and error messages only.
3	Show error messages only. This is the least verbose level.

You can call `libct_set_log_level()` from the application code to set the log level.

Log messages are written to standard output (stdout) and will be printed to the terminal window where the application was started by default. So to view the log messages, the application must be started from the terminal window explicitly. The logs can be viewed and saved as follows assuming *appExe* is the application executable.

```
appExe | tee libcaretaker.log
```


Chapter 3

Monitoring the Device

This section documents monitoring Caretaker data after integrating the application with the library.

Overview

Two groups of APIs are defined to simplify getting started: Primary and Secondary APIs. The primary API is the core interface required to connect to the Caretaker device to monitor numeric and waveform data, and the secondary API is an auxiliary interface to parse, read and write additional Caretaker data.

The sequence diagram below illustrates how a user application will interact with the primary functions to connect and start monitoring Caretaker data, which can be summarized as the following six steps.

- **Step 1:** Initialize a library context to associate with a Caretaker device.
- **Step 2:** Discover the device.
- **Step 3:** Connect to the device.
- **Step 4:** Start monitoring device data.
- **Step 5:** Calibrate and start measurements.
- **Step 6:** Handle numeric and waveform data updates.




Figure 3.1 Sequence diagram to connect and monitor data.

Note: For simplification the code snippets used throughout this manual omit error handling. It is recommended that any application handle errors returned from library functions appropriately.

Step 1: Initialize a library context

Start by initializing a library context (or library instance) by calling `libct_init()`. Specify the class of device to associate with this library context in the init data. Set the appropriate application callback functions, and set any unused callbacks to null. At a minimum, the following callback functions should be implemented.

- `on_device_discovered()`
- `on_device_connected_ready()`
- `on_device_disconnected()`
- `on_data_received()`

Next call `libct_init()` with the context pointer, init data, and callback variables. Note the context pointer must be initialized to null prior to being passed to `libct_init()` to indicate it is not in use, otherwise `libct_init()` will return error.

```
libct_init_data_t init_data;
memset(&init_data, 0, sizeof(init_data));
init_data.device_class = LIBCT_DEVICE_CLASS_BLE_CARETAKER4;

libct_app_callbacks_t callbacks = { ... }

libct_context_t* context = NULL;
int status = libct_init(&context, &init_data, &callbacks);
if ( LIBCT_FAILED(status) ) {
    // Handle error
}
```

Optionally, after initializing the context, the application specific data can be saved in the context for retrieval and use later in the application callbacks. For example, a C++ main application instance can be set as app specific data for access in the callbacks. See `libct_set_app_specific_data()` and `libct_get_app_specific_data()` for more information.

```
libct_set_app_specific_data(context, this);
```

Step 2: Discover a device

If `libct_init()` returned success, a device context has been created and initialized to start device discovery. Call `libct_start_discovery()` to scan for nearby devices. It will scan for some specified timeout (20 seconds in the code example below) and automatically stop if the time out has been reached or if `libct_stop_discovery()` has been explicitly called to stop device discovery.

```
libct_start_discovery(context, 20000);
```

The application will receive notifications later from devices matching the device class specified in the `init_data` passed to `libct_init()`. These notifications will be signalled to the application with the following callback functions.

- `on_device_discovered()`
- `on_discovery_timedout()`
- `on_discovery_failed()`

Note: The `on_device_discovered()` callback must be implemented to receive notification when a matching device is found.

Step 3: Connect to the device

The following implementation illustrates connecting to the first device found. However, by implementing a device whitelist to check discovered devices against a known acceptable list, a specific device can be searched for and automatically connected or all discovered devices can be displayed on the application GUI allowing the user to select the appropriate device.

```
void LIBCTAPI on_device_discovered_cb(libct_context_t* context,
    libct_device_t* device) {
    libct_stop_discovery(context);
    libct_connect(context, device);
}
```

After calling `libct_connect()`, the application will later receive notifications signalling the connection status with the following callbacks.

- `on_device_connected_not_ready()`
- `on_device_connected_ready()`
- `on_connect_error()`
- `on_connect_timedout()`

Note: The `on_device_connected_ready()` callback must be implemented to receive notification when the connection is established and the device is ready to receive requests.

Step 4: Monitor device data

After the connection has been established, monitoring vitals from the device can be started. The following is a sample implementation illustrating this in the connection ready callback, but note monitoring the device data can be deferred until other application events are received, such as user input from the application GUI.

```
void LIBCTAPI on_device_connected_ready_cb(libct_context_t* context,
    libct_device_t* device) {
    int flags = (LIBCT_MONITOR_INT_PULSE |
        LIBCT_MONITOR_PARAM_PULSE |
        LIBCT_MONITOR_VITALS |
        LIBCT_MONITOR_CUFF_PRESSURE |
        LIBCT_MONITOR_DEVICE_STATUS |
        LIBCT_MONITOR_BATTERY_INFO);

    libct_start_monitoring(context, flags);
}
```

Note: The application will receive monitor status via the `on_start_monitoring()` callback. The callback will be invoked only once in response to each `libct_start_monitoring()` call and is thus a oneshot callback.

Step 5: Calibrate and start measurements

After calling `libct_start_monitoring()`, the application will start receiving data from the device via the `on_data_received()` callback, however, the application will not receive valid vitals and waveform data until the blood pressure measurements are calibrated. Again, starting calibration can be deferred until other application events are received, such as user input from the application GUI.

The following code illustrates starting automatic calibration. Note the patient posture must be retrieved elsewhere, such as from the application GUI.

```
libct_cal_t cal;
cal.type = LIBCT_AUTO_CAL;
cal.config.auto_cal.posture = posture;
libct_start_measuring(context, &cal);
```

And the following code illustrates starting manual calibration. Again, the systolic and diastolic initial values must be retrieved elsewhere, such as from the application GUI.

```
libct_cal_t cal;
cal.type = LIBCT_MANUAL_CAL;
cal.config.manual_cal.settings.systolic = systolic;
cal.config.manual_cal.settings.diastolic = diastolic;
libct_start_measuring(context, &cal);
```

Note: The application will receive measurement status via the `on_start_measuring()` callback, which is a oneshot callback, i.e., the callback will be invoked only once in response to each `libct_start_measuring()` call.

Step 6: Handle numeric and waveform data updates

If monitoring and measurements were started successfully, the application will start receiving numeric and waveform data updates. The application `on_data_received()` callback will be notified continuously while data is received from the device.

The following code snippet illustrates processing data received from the device in the application `on_data_received()` callback. See the `stream data` structure for data format details.

```
void LIBCTAPI on_data_received_cb(libct_context_t* context,
    libct_device_t* device, libct_stream_data_t* data) {
    // Obtain the application instance set earlier with libct_set_app_specific_data().
    // Note libct_get_app_specific_data() returns null if libct_set_app_specific_data() was not
    // called earlier to set the application instance.
    MainWindow* window = (MainWindow*) libct_get_app_specific_data(context);

    // Update device status
    if ( data->device_status->valid ) {
        // ... check device status flags
    }

    // Update vitals
    libct_vitals_t* vitals = libct_get_last_dp(data, vitals);
    if ( vitals && vitals->valid ) {
        window->setHr(vitals->heart_rate);
        window->setRes(vitals->respiration);
        window->setMap(vitals->map);
        window->setBp(vitals->systolic, vitals->diastolic);
    }

    // Update the pulse rate waveform
    unsigned int idx;
    libct_pulse_t* pulse;
    for_each_dp(data, idx, pulse, raw_pulse) {
        if ( pulse && pulse->valid ) {
            window->rawPulseWaveform->add(pulse->timestamp, pulse->
                value);
        }
    }
}
```

```
    }

    // Update the pulse pressure waveform
    for_each_dp(data, idx, pulse, int_pulse) {
        if ( pulse && pulse->valid ) {
            if ( pulse && pulse->valid ) {
                window->intPulseWaveform->add(pulse->timestamp, pulse->
                    value);
            }
        }
    }
}
```

Chapter 4

Module Index

4.1 Modules

Here is a list of all modules:

Device Information	17
Primary API	23
Secondary API	29

Chapter 5

Data Structure Index

5.1 Data Structures

Here are the data structures with brief descriptions:

libct_app_callbacks_t	Structure used to provide asynchronous notifications to the application	45
libct_battery_info_t	Battery info data point within the libct_stream_data_t packet	66
libct_bp_settings_t	Structure to write the caretaker manual blood pressure settings	67
libct_cal_curve_t	Calibration curve data point within the libct_stream_data_t packet	68
libct_cal_t	Structure used to pass calibration data to libct_start_measuring()	69
libct_cal_type_t	The Caretaker calibration types	71
libct_context_t	An opaque type representing a library instance associated with (or bound to) a device the application is monitoring	71
libct_cuff_pressure_t	Cuff pressure data point within the libct_stream_data_t packet	71
libct_device_class_t	Classes of devices that can be monitored by this library	73
libct_device_state_t	The Caretaker device states	73
libct_device_status_t	Device status data point within the libct_stream_data_t packet	74
libct_device_t	Handle used to identify a connected device the application is monitoring	78
libct_init_data_t	Structure defining initialization data passed to libct_init()	84
libct_monitor_flags_t	Data monitor flags passed to libct_start_monitoring()	85
libct_param_pulse_t	Parametrized pulse data within the libct_stream_data_t packet	85
libct_posture_t	Patient postures	89
libct_pulse_ox_t	Pulse oximetry data point within the libct_stream_data_t packet	89

libct_pulse_t	
Raw or integrated pulse data point within the libct_stream_data_t packet	90
libct_pulse_waveform_t	
Pulse waveform data returned from the device as a result of a previous read request	91
libct_status_t	
Function return status codes	92
libct_stream_data_t	
This structure is used to hand-off data received from the remote device to the application . . .	92
libct_temperature_t	
Temperature data point within the libct_stream_data_t packet	101
libct_version_t	
CareTaker version information	102
libct_vitals2_t	
Secondary Vitals data point within the libct_stream_data_t packet	102
libct_vitals_t	
Vitals data point within the libct_stream_data_t packet	104

Chapter 6

Module Documentation

6.1 Device Information

This module describes the interface to retrieve general information about the Caretaker device.

Data Structures

- class `libct_device_t`
Handle used to identify a connected device the application is monitoring.

Macros

- `#define libct_device_get_state(dev) (dev)->get_state(dev)`
Convenience macro to `device->get_state(device)`.
- `#define libct_device_uninitialized(dev) (((dev)->get_state(dev)) & LIBCT_STATE_UNINITIALIZED)`
Returns non-zero (true) if the device is not initialized, and zero (false) otherwise.
- `#define libct_device_initialized(dev) (!libct_device_uninitialized(dev))`
Returns non-zero (true) if the device is initialized, and zero (false) otherwise.
- `#define libct_device_discovering(dev) (((dev)->get_state(dev)) & LIBCT_STATE_DISCOVERING)`
Returns non-zero (true) if discovering the device, and zero (false) otherwise.
- `#define libct_device_connecting(dev) (((dev)->get_state(dev)) & LIBCT_STATE_CONNECTING)`
Returns non-zero (true) if connecting to the device, and zero (false) otherwise.
- `#define libct_device_connected(dev) (((dev)->get_state(dev)) & LIBCT_STATE_CONNECTED)`
Returns non-zero (true) if connected to the device, and zero (false) otherwise.
- `#define libct_device_disconnecting(dev) (((dev)->get_state(dev)) & LIBCT_STATE_DISCONNECTING)`
Returns non-zero (true) if disconnecting from the device, and zero (false) otherwise.
- `#define libct_device_disconnected(dev) (((dev)->get_state(dev)) & LIBCT_STATE_DISCONNECTED)`
Returns non-zero (true) if disconnected from the device, and zero (false) otherwise.
- `#define libct_device_monitoring(dev) (((dev)->get_state(dev)) & LIBCT_STATE_MONITORING)`
Returns non-zero (true) if receiving data from the device, and zero (false) otherwise.
- `#define libct_device_measuring(dev) (((dev)->get_state(dev)) & LIBCT_STATE_MEASURING)`
Returns non-zero (true) if taking blood pressure measurements, and zero (false) otherwise.
- `#define libct_device_get_class(dev) (dev)->get_class(dev)`
Convenience macro to `device->get_class(device)`.

- `#define libct_device_get_name(dev) (dev)->get_name(dev)`
Convenience macro to `device->get_name(device)`.
- `#define libct_device_get_address(dev) (dev)->get_address(dev)`
Convenience macro to `device->get_address(device)`.
- `#define libct_device_get_serial_number(dev) (dev)->get_serial_number(dev)`
Convenience macro to `device->get_serial_number(device)`.
- `#define libct_device_get_hw_version(dev) (dev)->get_hw_version(dev)`
Convenience macro to `device->get_hw_version(device)`.
- `#define libct_device_get_fw_version(dev) (dev)->get_fw_version(dev)`
Convenience macro to `device->get_fw_version(device)`.
- `#define libct_device_get_context(dev) (dev)->get_context(dev)`
Convenience macro to `device->get_context(device)`.

6.1.1 Detailed Description

This module describes the interface to retrieve general information about the Caretaker device.

6.1.2 Macro Definition Documentation

6.1.2.1 libct_device_get_state

```
#define libct_device_get_state(  
    dev ) (dev)->get_state(dev)
```

Convenience macro to `device->get_state(device)`.

Parameters

<code>dev</code>	Pointer to the device instance.
------------------	---------------------------------

6.1.2.2 libct_device_uninitialized

```
#define libct_device_uninitialized(  
    dev ) (((dev)->get_state(dev)) & LIBCT_STATE_UNINITIALIZED)
```

Returns non-zero (true) if the device is not initialized, and zero (false) otherwise.

Parameters

<code>dev</code>	Pointer to device instance.
------------------	-----------------------------

6.1.2.3 libct_device_discovering

```
#define libct_device_discovering(  
    dev ) (((dev)->get_state(dev)) & LIBCT_STATE_DISCOVERING)
```

Returns non-zero (true) if discovering the device, and zero (false) otherwise.

Parameters

<i>dev</i>	Pointer to device instance.
------------	-----------------------------

6.1.2.4 libct_device_connecting

```
#define libct_device_connecting(  
    dev ) (((dev)->get_state(dev)) & LIBCT_STATE_CONNECTING)
```

Returns non-zero (true) if connecting to the device, and zero (false) otherwise.

Parameters

<i>dev</i>	Pointer to device instance.
------------	-----------------------------

6.1.2.5 libct_device_connected

```
#define libct_device_connected(  
    dev ) (((dev)->get_state(dev)) & LIBCT_STATE_CONNECTED)
```

Returns non-zero (true) if connected to the device, and zero (false) otherwise.

Parameters

<i>dev</i>	Pointer to device instance.
------------	-----------------------------

6.1.2.6 libct_device_disconnecting

```
#define libct_device_disconnecting(  
    dev ) (((dev)->get_state(dev)) & LIBCT_STATE_DISCONNECTING)
```

Returns non-zero (true) if disconnecting from the device, and zero (false) otherwise.

Parameters

<i>dev</i>	Pointer to device instance.
------------	-----------------------------

6.1.2.7 libct_device_disconnected

```
#define libct_device_disconnected(  
    dev ) (((dev)->get_state(dev)) & LIBCT_STATE_DISCONNECTED)
```

Returns non-zero (true) if disconnected from the device, and zero (false) otherwise.

Parameters

<i>dev</i>	Pointer to device instance.
------------	-----------------------------

6.1.2.8 libct_device_monitoring

```
#define libct_device_monitoring(  
    dev ) (((dev)->get_state(dev)) & LIBCT_STATE_MONITORING)
```

Returns non-zero (true) if receiving data from the device, and zero (false) otherwise.

Parameters

<i>dev</i>	Pointer to device instance.
------------	-----------------------------

6.1.2.9 libct_device_measuring

```
#define libct_device_measuring(  
    dev ) (((dev)->get_state(dev)) & LIBCT_STATE_MEASURING)
```

Returns non-zero (true) if taking blood pressure measurements, and zero (false) otherwise.

Parameters

<i>dev</i>	Pointer to device instance.
------------	-----------------------------

6.1.2.10 libct_device_get_class

```
#define libct_device_get_class(  
    dev ) (dev)->get_class(dev)
```

Convenience macro to [device->get_class\(device\)](#).

Parameters

<i>dev</i>	Pointer to the device instance.
------------	---------------------------------

6.1.2.11 libct_device_get_name

```
#define libct_device_get_name(  
    dev ) (dev)->get_name(dev)
```

Convenience macro to [device->get_name\(device\)](#).

Parameters

<i>dev</i>	Pointer to the device instance.
------------	---------------------------------

6.1.2.12 libct_device_get_address

```
#define libct_device_get_address(  
    dev ) (dev)->get_address(dev)
```

Convenience macro to [device->get_address\(device\)](#).

Parameters

<i>dev</i>	Pointer to the device instance.
------------	---------------------------------

6.1.2.13 libct_device_get_serial_number

```
#define libct_device_get_serial_number(  
    dev ) (dev)->get_serial_number(dev)
```

Convenience macro to [device->get_serial_number\(device\)](#).

Parameters

<i>dev</i>	Pointer to the device instance.
------------	---------------------------------

6.1.2.14 libct_device_get_hw_version

```
#define libct_device_get_hw_version(  
    dev ) (dev)->get_hw_version(dev)
```

Convenience macro to [device->get_hw_version\(device\)](#).

Parameters

<i>dev</i>	Pointer to the device instance.
------------	---------------------------------

6.1.2.15 libct_device_get_fw_version

```
#define libct_device_get_fw_version(  
    dev ) (dev)->get_fw_version(dev)
```

Convenience macro to [device->get_fw_version\(device\)](#).

Parameters

<i>dev</i>	Pointer to the device instance.
------------	---------------------------------

6.1.2.16 libct_device_get_context

```
#define libct_device_get_context(  
    dev ) (dev)->get_context(dev)
```

Convenience macro to [device->get_context\(device\)](#).

Parameters

<i>dev</i>	Pointer to the device instance.
------------	---------------------------------

6.2 Primary API

The group of primary functions that are required to connect to a Caretaker device and monitor data.

Functions

- LIBCTEXTEXPORT int `libct_init` (`libct_context_t` **context, `libct_init_data_t` *data, `libct_app_callbacks_t` *callbacks)
Initializes device context.
- LIBCTEXTEXPORT void `libct_deinit` (`libct_context_t` *context)
De-initializes the context.
- LIBCTEXTEXPORT int `libct_start_discovery` (`libct_context_t` *context, unsigned long timeout)
Discover the device.
- LIBCTEXTEXPORT int `libct_stop_discovery` (`libct_context_t` *context)
Stop device discovery.
- LIBCTEXTEXPORT int `libct_connect` (`libct_context_t` *context, `libct_device_t` *device)
Connect to a device.
- LIBCTEXTEXPORT int `libct_disconnect` (`libct_context_t` *context)
Disconnect from a device.
- LIBCTEXTEXPORT int `libct_start_monitoring` (`libct_context_t` *context, int flags)
Start monitoring data at the remote caretaker device.
- LIBCTEXTEXPORT int `libct_stop_monitoring` (`libct_context_t` *context)
Stops monitoring.
- LIBCTEXTEXPORT int `libct_start_measuring` (`libct_context_t` *context, `libct_cal_t` *cal)
Start taking measurement.
- LIBCTEXTEXPORT int `libct_stop_measuring` (`libct_context_t` *context)
Stops measuring.

6.2.1 Detailed Description

The group of primary functions that are required to connect to a Caretaker device and monitor data.

6.2.2 Function Documentation

6.2.2.1 `libct_init()`

```
LIBCTEXTEXPORT int libct_init (
    libct_context_t ** context,
    libct_init_data_t * data,
    libct_app_callbacks_t * callbacks )
```

Initializes device context.

Call this function to initialize a device context before calling any other library functions with the said context.

Note

You can initialize multiple contexts if you wish to connect to multiple devices simultaneously, but you must call `libct_deinit()` from the same thread to de-initialize each context when it is no longer needed.

Parameters

<i>context</i>	Address to store the created context. IMPORTANT: Initialize the context pointer to null before passing it. The internal library code depends on this to ensure the context is initialized only once.
<i>data</i>	Data to initialize the context.
<i>callbacks</i>	The application callback functions to receive asynchronous notifications. This pointer must not be null, or else your application will not receive notifications notifying connection and data events. However, you can set function pointers within this structure that you don't care about to null. NOTE: You can set application specific data to use inside your callbacks with libct_set_app_specific_data() after initialization, and later retrieve it with libct_get_app_specific_data() to get the application instance data to act upon inside the callbacks.

Returns

An appropriate [status code](#) indicating success or error.

6.2.2.2 libct_deinit()

```
LIBCTEXPORT void libct_deinit (
    libct_context_t * context )
```

De-initializes the context.

Call this function to release resources when you no longer need the context.

IMPORTANT: The application must call `libct_deinit()` some time after calling [libct_init\(\)](#) to prevent resource leaks. [libct_deinit\(\)](#) must not be called from any library callback function. Library callbacks are called from internal library threads that this function attempts to kill. As such, it must only be called from an application thread.

Parameters

<i>context</i>	The context returned from libct_init() .
----------------	--

6.2.2.3 libct_start_discovery()

```
LIBCTEXPORT int libct_start_discovery (
    libct_context_t * context,
    unsigned long timeout )
```

Discover the device.

This function start scanning for devices specify by the [device class](#) in the initialization data passed earlier to [libct_init\(\)](#). Scan results will be notified asynchronously via the [application callbacks](#) passed to [libct_init\(\)](#); specifically, these discovery callback functions will be invoked some time later with the results when devices are discovered or if scanning timed out or failed.

- [on_device_discovered\(\)](#)
- [on_discovery_timeout\(\)](#)
- [on_discovery_failed\(\)](#)

Note

Devices must be advertising and be within range for this method to discover them. Press the button on the caretaker to start advertising. Note the caretaker only advertises for 20 seconds after pressing the button and then stops.

Parameters

<i>context</i>	The context returned from libct_init() .
<i>timeout</i>	Scanning will be canceled after the number of milliseconds specified by this timeout and the application discovery timeout callback will be invoked.

Returns

An appropriate [status code](#) indicating success or error.

6.2.2.4 libct_stop_discovery()

```
LIBCTEXPORT int libct_stop_discovery (
    libct_context_t * context )
```

Stop device discovery.

Call this function to stop device discovery previously started with [libct_start_discovery\(\)](#).

Parameters

<i>context</i>	The context returned from libct_init() .
----------------	--

Returns

An appropriate [status code](#) indicating success or error.

6.2.2.5 libct_connect()

```
LIBCTEXPORT int libct_connect (
    libct_context_t * context,
    libct_device_t * device )
```

Connect to a device.

Call this method after device discovery to establish connection to the device. The results will be notified asynchronously via the [application callbacks](#) passed to [libct_init\(\)](#); specifically, one or more of the following callback functions will be invoked some time later with the results if the connection is established, timed out or failed.

- [on_device_connected_not_ready\(\)](#)
- [on_device_connected_ready\(\)](#)
- [on_connect_error\(\)](#)
- [on_connect_timedout\(\)](#)

Parameters

<i>context</i>	The context returned from libct_init() .
<i>device</i>	A discovered device notified with your application callback on_device_discovered() .

Returns

An appropriate [status code](#) indicating success or error.

6.2.2.6 libct_disconnect()

```
LIBCTEXPORT int libct_disconnect (
    libct_context_t * context )
```

Disconnect from a device.

Call this method after calling [libct_connect\(\)](#) to clean up resources that were allocated by the connect call.

IMPORTANT: The application must call [libct_disconnect\(\)](#) to release connection resources before calling [libct_connect\(\)](#) subsequently on the same context. Otherwise, the subsequent connect calls may fail. Also, the application must not call [libct_disconnect\(\)](#) from any library callback function. Library callbacks are called from internal library threads and this function attempts to kill. As such, it must only be called from an application thread.

Parameters

<i>context</i>	The context returned from libct_init() .
----------------	--

Returns

An appropriate [status code](#) indicating success or error.

6.2.2.7 libct_start_monitoring()

```
LIBCTEXPORT int libct_start_monitoring (
    libct_context_t * context,
    int flags )
```

Start monitoring data at the remote caretaker device.

Call this function after the connection is established with the device to start monitoring data or to change the data being monitored.

Calling this function will trigger your application's [on_start_monitoring\(\)](#) to be invoked some time later with results. Also, if monitoring was started successfully, data from the device will be notified to your application's [on_data_received\(\)](#) continuously until stopped explicitly by calling to [libct_stop_monitoring\(\)](#) or if the device becomes disconnected.

Parameters

<i>context</i>	The context returned from libct_init() .
<i>flags</i>	Bitwise OR of monitor flags specifying the data to monitor. Note: The stream data packets notified to the application depends on these flags. So you can control the amount of data reported to the application by specifying only the monitoring flags corresponding to the data you care about.

Returns

An appropriate [status code](#) indicating success or error.

6.2.2.8 libct_stop_monitoring()

```
LIBCTEXPORT int libct_stop_monitoring (
    libct_context_t * context )
```

Stops monitoring.

Call this method to stop monitoring data after calling [libct_start_monitoring\(\)](#) successfully.

Calling this function will trigger [on_stop_monitoring\(\)](#) to be invoked sometime later with results.

Parameters

<i>context</i>	The context returned from libct_init() .
----------------	--

Returns

An appropriate [status code](#) indicating success or error.

6.2.2.9 libct_start_measuring()

```
LIBCTEXPORT int libct_start_measuring (
    libct_context_t * context,
    libct_cal_t * cal )
```

Start taking measurement.

If monitoring was started successfully, call this function to initialize (calibrate) blood pressure settings with either auto or manual calibration then start taking vital sign measurements.

Calling this function will trigger [on_start_measuring\(\)](#) to be invoked sometime later with results.

Parameters

<i>context</i>	The context returned from libct_init() .
<i>cal</i>	Auto or manual calibration settings.

Returns

An appropriate [status code](#) indicating success or error

6.2.2.10 libct_stop_measuring()

```
LIBCTEXPORT int libct_stop_measuring (
    libct_context_t * context )
```

Stops measuring.

Call this method to stop measuring data after calling [libct_start_measuring\(\)](#) successfully.

Calling this function will trigger [on_stop_measuring\(\)](#) to be invoked sometime later with results.

Parameters

<i>context</i>	The context returned from libct_init() .
----------------	--

Returns

An appropriate [status code](#) indicating success or error.

6.3 Secondary API

The group of auxiliary functions and macros available to parse, read and write additional Caretaker device data.

Macros

- `#define libct_dp_count(data, memb) (data)->memb.count`
Returns the count of data points of the specified member array contained in [stream data](#) received at the application.
- `#define libct_get_dp(data, memb, pos)`
Extract a single data point from the specified member array contained in [stream data](#) received at the application.
- `#define libct_get_last_dp(data, memb) libct_get_dp(data, memb, (data)->memb.count-1)`
Extract the newest data point from the specified member array contained in [stream data](#) received at the application.
- `#define libct_get_first_dp(data, memb) libct_get_dp(data, memb, 0)`
Extract the oldest data point from the specified member array contained in [stream data](#) received at the application.
- `#define for_each_dp(data, idx, dp, memb) for(idx=0; (idx<(data)->memb.count) && (dp=((data)->memb.↵datapoints) ? &(data)->memb.datapoints[idx] : NULL); idx++)`
Iterate over data points of the specified member array to extract from [stream data](#) received at the application.
- `#define libct_inc_cuff_pressure(context) libct_adjust_cuff_pressure((context), 1)`
Increments the cuff pressure in 10 mmHg increment.
- `#define libct_dec_cuff_pressure(context) libct_adjust_cuff_pressure((context), 0)`
Decrements the cuff pressure in 10 mmHg increment.

Functions

- `LIBCTEXPORT libct_device_t * libct_get_device (libct_context_t *context)`
Returns the device handle.
- `LIBCTEXPORT void libct_set_app_specific_data (libct_context_t *context, void *data)`
Sets application specific data that can be retrieved and used later in the callbacks.
- `LIBCTEXPORT void * libct_get_app_specific_data (libct_context_t *context)`
Retrieve application specific data.
- `LIBCTEXPORT const char * libct_get_version_string (void)`
Get the library version info.
- `LIBCTEXPORT const char * libct_get_build_date_string (void)`
Get the library build date and time string.
- `LIBCTEXPORT void libct_set_log_level (int level)`
Sets the library log level.
- `LIBCTEXPORT int libct_recalibrate (libct_context_t *context)`
Re-calibrates the device.
- `LIBCTEXPORT int libct_adjust_cuff_pressure (libct_context_t *context, int direction)`
Adjusts the cuff pressure in 10 mmHg increment/decrement.
- `LIBCTEXPORT int libct_rd_cuff_pressure (libct_context_t *context)`
Reads the cuff pressure.
- `LIBCTEXPORT int libct_vent_cuff (libct_context_t *context)`
Deflates the cuff pressure.
- `LIBCTEXPORT int libct_clr_status (libct_context_t *context)`
Clears the device status.
- `LIBCTEXPORT int libct_diag_flush (libct_context_t *context)`
Invoke the device diagnostic plumbing tree flush.
- `LIBCTEXPORT int libct_wrt_snr_min (libct_context_t *context, int snr)`

- Writes the device noise filter parameter.*
- LIBCTEXPORT int [libct_rd_snr_min](#) ([libct_context_t](#) *context)
- Reads the device noise filter parameter.*
- LIBCTEXPORT int [libct_wrt_display_state](#) ([libct_context_t](#) *context, unsigned char state)
- Turns the device display on/off.*
- LIBCTEXPORT int [libct_rd_display_state](#) ([libct_context_t](#) *context)
- Reads the device display state.*
- LIBCTEXPORT int [libct_wrt_recal_itvl](#) ([libct_context_t](#) *context, unsigned int itvl)
- Writes the recalibration interval.*
- LIBCTEXPORT int [libct_rd_recal_itvl](#) ([libct_context_t](#) *context)
- Reads the recalibration interval.*
- LIBCTEXPORT int [libct_wrt_waveform_clamping](#) ([libct_context_t](#) *context, unsigned char value)
- Writes the device waveform clamping setting.*
- LIBCTEXPORT int [libct_rd_waveform_clamping](#) ([libct_context_t](#) *context)
- Reads the device waveform clamping setting.*
- LIBCTEXPORT int [libct_rd_vitals_filter](#) ([libct_context_t](#) *context)
- Reads the current filter settings to enable or disable filtering outlier vitals measurements.*
- LIBCTEXPORT int [libct_wrt_vitals_filter](#) ([libct_context_t](#) *context, unsigned char value)
- Writes the filter settings to enable or disable filtering outlier vitals measurements.*
- LIBCTEXPORT int [libct_wrt_simulation_mode](#) ([libct_context_t](#) *context, unsigned char mode)
- Writes the device simulation mode.*
- LIBCTEXPORT int [libct_wrt_motion_timeout](#) ([libct_context_t](#) *context, int timeout)
- Writes the motion tolerance timeout parameter.*
- LIBCTEXPORT int [libct_rd_motion_timeout](#) ([libct_context_t](#) *context)
- Reads the motion tolerance timeout parameter.*
- LIBCTEXPORT int [libct_rd_persistent_log](#) ([libct_context_t](#) *context)
- Reads the device log messages.*

6.3.1 Detailed Description

The group of auxiliary functions and macros available to parse, read and write additional Caretaker device data.

6.3.2 Macro Definition Documentation

6.3.2.1 libct_dp_count

```
#define libct_dp_count(  
    data,  
    memb ) (data)->memb.count
```

Returns the count of data points of the specified member array contained in [stream data](#) received at the application.

Parameters

<i>data</i>	The stream data received in your on_data_received() application callback.
<i>memb</i>	The stream data member whose data point count is being queried.

Returns

The extracted data point on success, and null on failure.

6.3.2.2 libct_get_dp

```
#define libct_get_dp(  
    data,  
    memb,  
    pos )
```

Value:

```
((\  
    __typeof__((data)->memb.datapoints[0]) *dp = NULL; \  
    if ( (data)->memb.count && (pos) < (data)->memb.count ) { \  
        dp = &(data)->memb.datapoints[(pos)]; \  
    } \  
    (dp); \  
))
```

Extract a single data point from the specified member array contained in [stream data](#) received at the application.

Parameters

<i>data</i>	The stream data received in your on_data_received() application callback.
<i>memb</i>	The stream data member name of the data point to extract.
<i>pos</i>	The position of the data point to extract.

Returns

The extracted data point on success, and null on failure.

6.3.2.3 libct_get_last_dp

```
#define libct_get_last_dp(  
    data,  
    memb ) libct_get_dp(data, memb, (data)->memb.count-1)
```

Extract the newest data point from the specified member array contained in [stream data](#) received at the application.

Parameters

<i>data</i>	The stream data received in your on_data_received() application callback.
<i>memb</i>	The stream data member name of the data point to extract.

Returns

The extracted data point on success, and null on failure.

6.3.2.4 libct_get_first_dp

```
#define libct_get_first_dp(  
    data,  
    memb ) libct_get_dp(data, memb, 0)
```

Extract the oldest data point from the specified member array contained in [stream data](#) received at the application.

Parameters

<i>data</i>	The stream data received in your on_data_received() application callback.
<i>memb</i>	The stream data member name of the data point to extract.

Returns

The extracted data point on success, and null on failure.

6.3.2.5 for_each_dp

```
#define for_each_dp(  
    data,  
    idx,  
    dp,  
    memb ) for(idx=0; (idx<(data)->memb.count) && (dp=((data)->memb.datapoints) ?  
&(data)->memb.datapoints[idx] : NULL); idx++)
```

Iterate over data points of the specified member array to extract from [stream data](#) received at the application.

Parameters

<i>data</i>	The stream data received in your on_data_received() application callback.
<i>idx</i>	Iterator variable of type unsigned integer.
<i>dp</i>	Pointer variable of type corresponding to the memb argument.
<i>memb</i>	The stream data member name of the data points to extract.

6.3.2.6 libct_inc_cuff_pressure

```
#define libct_inc_cuff_pressure(  
    context ) libct_adjust_cuff_pressure((context), 1)
```

Increments the cuff pressure in 10 mmHg increment.

Parameters

<i>context</i>	The context returned from libct_init() .
----------------	--

Returns

Success if a request was queued to be sent to the device, and error otherwise.

6.3.2.7 libct_dec_cuff_pressure

```
#define libct_dec_cuff_pressure(  
    context ) libct_adjust_cuff_pressure((context), 0)
```

Decrements the cuff pressure in 10 mmHg increment.

Parameters

<i>context</i>	The context returned from libct_init() .
----------------	--

Returns

Success if a request was queued to be sent to the device, and error otherwise.

6.3.3 Function Documentation

6.3.3.1 libct_get_device()

```
LIBCTEXPORT libct_device_t* libct_get_device (  
    libct_context_t * context )
```

Returns the device handle.

Call this function to get a pointer to the device instance associated with the context.

Parameters

<i>context</i>	The context returned from libct_init() .
----------------	--

Returns

Device object or NULL;

See also

The device APIs

6.3.3.2 libct_set_app_specific_data()

```
LIBCTEXPORT void libct_set_app_specific_data (
    libct_context_t * context,
    void * data )
```

Sets application specific data that can be retrieved and used later in the callbacks.

Basically, this function provides the means to bind your callback application code with the library context. For example, set application instance data after initializing the library, and then retrieve the instance data using [libct_get_app_specific_data\(\)](#) inside the [callbacks](#).

```
// QT main window initialization
MainWindow::MainWindow(QWidget *parent) :
    QMainWindow(parent),
    ui(new Ui::MainWindow)
{
    ui->setupUi(this);
    setWindowTitle(QString("SampleApp"));

    // initialize ui (code not shown)

    // initialize the library
    int status = libct_init(&context, &init_data, &callbacks);
    if ( LIBCT_FAILED(status) ) {
        // handle error
    }

    // set main window instance to act upon in the callbacks
    libct_set_app_specific_data(context, this);
}

// libcaretaker callback - called to notify data from the device
void on_data_received_cb(libct_context_t* context, libct_device_t* device,
    libct_stream_data_t* data) {
    MainWindow* window = (MainWindow*) libct_get_app_specific_data(context);

    // display the most recent vitals
    libct_vitals_t* vitals = libct_get_last_dp(data, vitals);
    if ( vitals && vitals->valid ) {
        window->setHr(vitals->heart_rate);
        window->setRes(vitals->respiration);
        window->setMap(vitals->map);
        window->setBp(vitals->systolic, vitals->diastolic);
    }
}
```

Parameters

<i>context</i>	The context returned from libct_init() .
<i>data</i>	Generic pointer to the application specific data, or null to clear the existing pointer.

6.3.3.3 libct_get_app_specific_data()

```
LIBCTEXPORT void* libct_get_app_specific_data (
    libct_context_t * context )
```

Retrieve application specific data.

Retrieves application specific data last set with [libct_set_app_specific_data\(\)](#).

Parameters

<i>context</i>	The context returned from libct_init() .
----------------	--

6.3.3.4 libct_get_version_string()

```
LIBCTEXPORT const char* libct_get_version_string (
    void )
```

Get the library version info.

Call this function to get the library version.

Returns

The library version string.

6.3.3.5 libct_get_build_date_string()

```
LIBCTEXPORT const char* libct_get_build_date_string (
    void )
```

Get the library build date and time string.

Call this function to get the library build date and time.

Returns

The library build date string.

6.3.3.6 libct_set_log_level()

```
LIBCTEXPORT void libct_set_log_level (
    int level )
```

Sets the library log level.

Call this function to set the log level to increase or decrease log messages verbosity.

Parameters

<i>level</i>	One of the following log levels. <ul style="list-style-type: none">• 0 - shows all logs, most verbose• 1 - shows only info, warning, and error logs• 2 - shows only warning and error logs• 3 - shows only error logs, least verbose
--------------	---

6.3.3.7 libct_recalibrate()

```
LIBCTEXPORT int libct_recalibrate (  
    libct_context_t * context )
```

Re-calibrates the device.

Call this function sometime later after calling [libct_start_measuring\(\)](#) to force vital signs re-calibration at the device while taking measurements.

Parameters

<i>context</i>	The context returned from libct_init() .
----------------	--

Returns

Success if a request was queued to be sent to the device, and error otherwise.

6.3.3.8 libct_adjust_cuff_pressure()

```
LIBCTEXPORT int libct_adjust_cuff_pressure (  
    libct_context_t * context,  
    int direction )
```

Adjusts the cuff pressure in 10 mmHg increment/decrement.

Note

The macros [libct_inc_cuff_pressure\(\)](#) and [libct_dec_cuff_pressure\(\)](#) simplify this function so you should use them instead.

Parameters

<i>context</i>	The context returned from libct_init() .
<i>direction</i>	Zero - Decrement pressure. Nonzero - Increment pressures

See also

[libct_inc_cuff_pressure\(\)](#)
[libct_dec_cuff_pressure\(\)](#)

Returns

Success if a request was queued to be sent to the device, and error otherwise.

6.3.3.9 libct_rd_cuff_pressure()

```
LIBCTEXPORT int libct_rd_cuff_pressure (
    libct_context_t * context )
```

Reads the cuff pressure.

Results will be notified later with [on_rd_cuff_pressure_rsp\(\)](#).

Parameters

<i>context</i>	The context returned from libct_init() .
----------------	--

Returns

Success if a request was queued to be sent to the device, and error otherwise.

6.3.3.10 libct_vent_cuff()

```
LIBCTEXPORT int libct_vent_cuff (
    libct_context_t * context )
```

Deflates the cuff pressure.

Results will be notified later with [on_vent_cuff_rsp\(\)](#).

Parameters

<i>context</i>	The context returned from libct_init() .
----------------	--

Returns

Success if a request was queued to be sent to the device, and error otherwise.

6.3.3.11 libct_clr_status()

```
LIBCTEXPORT int libct_clr_status (
    libct_context_t * context )
```

Clears the device status.

Results will be notified later with [on_clr_status_rsp\(\)](#).

Parameters

<i>context</i>	The context returned from libct_init() .
----------------	--

Returns

Success if a request was queued to be sent to the device, and error otherwise.

6.3.3.12 libct_diag_flush()

```
LIBCTEXPORT int libct_diag_flush (
    libct_context_t * context )
```

Invoke the device diagnostic plumbing tree flush.

Results will be notified later with [on_diag_flush_rsp\(\)](#).

Parameters

<i>context</i>	The context returned from libct_init() .
----------------	--

Returns

Success if a request was queued to be sent to the device, and error otherwise.

6.3.3.13 libct_wrt_snr_min()

```
LIBCTEXPORT int libct_wrt_snr_min (
    libct_context_t * context,
    int snr )
```

Writes the device noise filter parameter.

Results will be notified later with [on_wrt_snr_min_rsp\(\)](#).

Parameters

<i>context</i>	The context returned from libct_init() .
<i>snr</i>	The minimum signal-to-noise value. Valid range [0, 100].

Returns

Success if a request was queued to be sent to the device, and error otherwise.

6.3.3.14 libct_rd_snr_min()

```
LIBCTEXPORT int libct_rd_snr_min (  
    libct_context_t * context )
```

Reads the device noise filter parameter.

Results will be notified later with [on_rd_snr_min_rsp\(\)](#).

Parameters

<i>context</i>	The context returned from libct_init() .
----------------	--

Returns

Success if a request was queued to be sent to the device, and error otherwise.

6.3.3.15 libct_wrt_display_state()

```
LIBCTEXPORT int libct_wrt_display_state (  
    libct_context_t * context,  
    unsigned char state )
```

Turns the device display on/off.

Results will be notified later with [on_rd_snr_min_rsp\(\)](#).

Parameters

<i>context</i>	The context returned from libct_init() .
<i>state</i>	Display state to write: 0 = off, 1 = on.

Returns

Success if a request was queued to be sent to the device, and error otherwise.

6.3.3.16 libct_rd_display_state()

```
LIBCTEXPORT int libct_rd_display_state (
    libct_context_t * context )
```

Reads the device display state.

Results will be notified later with [on_rd_snr_min_rsp\(\)](#).

Parameters

<i>context</i>	The context returned from libct_init() .
----------------	--

Returns

Success if a request was queued to be sent to the device, and error otherwise.

6.3.3.17 libct_wrt_recal_itvl()

```
LIBCTEXPORT int libct_wrt_recal_itvl (
    libct_context_t * context,
    unsigned int itvl )
```

Writes the recalibration interval.

Results will be notified later with [on_wrt_recal_itvl_rsp\(\)](#).

Parameters

<i>context</i>	The context returned from libct_init() .
<i>itvl</i>	The recalibration interval in minutes. The acceptable range is [30, 240].

Returns

Success if a request was queued to be sent to the device, and error otherwise.

6.3.3.18 libct_rd_recal_itvl()

```
LIBCTEXPORT int libct_rd_recal_itvl (
    libct_context_t * context )
```

Reads the recalibration interval.

Results will be notified later with [on_rd_recal_itvl_rsp\(\)](#).

Parameters

<i>context</i>	The context returned from libct_init() .
----------------	--

Returns

Success if a request was queued to be sent to the device, and error otherwise.

6.3.3.19 libct_wrt_waveform_clamping()

```
LIBCTEXPORT int libct_wrt_waveform_clamping (
    libct_context_t * context,
    unsigned char value )
```

Writes the device waveform clamping setting.

Status will be notified later with [on_wrt_waveform_clamping\(\)](#).

Parameters

<i>context</i>	The context returned from libct_init() .
<i>value</i>	Clamp setting: 1 = ON, 0 = OFF.

Returns

Success if a request was queued to be sent to the device, and error otherwise.

6.3.3.20 libct_rd_waveform_clamping()

```
LIBCTEXPORT int libct_rd_waveform_clamping (
    libct_context_t * context )
```

Reads the device waveform clamping setting.

Status will be notified later with [on_rd_waveform_clamping\(\)](#).

Parameters

<i>context</i>	The context returned from libct_init() .
----------------	--

Returns

Success if a request was queued to be sent to the device, and error otherwise.

6.3.3.21 libct_rd_vitals_filter()

```
LIBCTEXPORT int libct_rd_vitals_filter (
    libct_context_t * context )
```

Reads the current filter settings to enable or disable filtering outlier vitals measurements.

Status will be notified later with [on_rd_vitals_filter\(\)](#).

Parameters

<i>context</i>	The context returned from libct_init() .
----------------	--

Returns

Success if a request was queued to be sent to the device, and error otherwise.

6.3.3.22 libct_wrt_vitals_filter()

```
LIBCTEXPORT int libct_wrt_vitals_filter (
    libct_context_t * context,
    unsigned char value )
```

Writes the filter settings to enable or disable filtering outlier vitals measurements.

Status will be notified later with [on_wrt_vitals_filter\(\)](#).

Parameters

<i>context</i>	The context returned from libct_init() .
<i>value</i>	Filter setting value: 1 = Enable, 0 = Disable.

Returns

Success if a request was queued to be sent to the device, and error otherwise.

6.3.3.23 libct_wrt_simulation_mode()

```
LIBCTEXPORT int libct_wrt_simulation_mode (
    libct_context_t * context,
    unsigned char mode )
```

Writes the device simulation mode.

Note

The device does not provide real-time data when simulation mode is enabled. Hard-coded numeric and waveform data (i.e., synthetic data) is provided. As such, simulation mode should be enabled for demonstration and test purposes only.

Parameters

<i>context</i>	The context returned from libct_init() .
<i>mode</i>	Simulation mode: 1 = Enable, 0 = Disable.

Returns

Success if a request was queued to be sent to the device, and error otherwise.

6.3.3.24 libct_wrt_motion_timeout()

```
LIBCTEXPORT int libct_wrt_motion_timeout (
    libct_context_t * context,
    int timeout )
```

Writes the motion tolerance timeout parameter.

Parameters

<i>context</i>	The context returned from libct_init() .
<i>timeout</i>	Time out in seconds. Acceptable range [0, 30]

Returns

Success if a request was queued to be sent to the device, and error otherwise.

6.3.3.25 libct_rd_motion_timeout()

```
LIBCTEXPORT int libct_rd_motion_timeout (
    libct_context_t * context )
```

Reads the motion tolerance timeout parameter.

Status will be notified later with [on_rd_motion_timeout\(\)](#).

Parameters

<i>context</i>	The context returned from libct_init() .
----------------	--

Returns

Success if a request was queued to be sent to the device, and error otherwise.

6.3.3.26 libct_rd_persistent_log()

```
LIBCTEXPORT int libct_rd_persistent_log (  
    libct_context_t * context )
```

Reads the device log messages.

Status will be notified later with [on_rd_persistent_log\(\)](#).

Note

Reading the device log is a slow request so the results will be delayed by many seconds.

Parameters

<i>context</i>	The context returned from libct_init() .
----------------	--

Returns

Success if a request was queued to be sent to the device, and error otherwise.

Chapter 7

Data Structure Documentation

7.1 libct_app_callbacks_t Struct Reference

Structure used to provide asynchronous notifications to the application.

Data Fields

- void(* [on_device_discovered](#))(libct_context_t *context, libct_device_t *device)
Function pointer to the application callback to receive scan notifications in response to calling [libct_start_discovery\(\)](#).
- void(* [on_discovery_timedout](#))(libct_context_t *context)
Function pointer to the application callback to receive timeout notification in response to calling [libct_start_discovery\(\)](#).
- void(* [on_discovery_failed](#))(libct_context_t *context, int error)
Function pointer to the application callback to receive error notification in response to calling [libct_start_discovery\(\)](#).
- void(* [on_device_connected_not_ready](#))(libct_context_t *context, libct_device_t *device)
Function pointer to the application callback to receive early connection notification in response to calling [libct_connect\(\)](#).
- void(* [on_device_connected_ready](#))(libct_context_t *context, libct_device_t *device)
Function pointer to the application callback to receive connection notification in response to calling [libct_connect\(\)](#).
- void(* [on_connect_error](#))(libct_context_t *context, libct_device_t *device, const char *error)
Function pointer to the application callback to receive error notification in response to calling [libct_connect\(\)](#).
- void(* [on_connect_timedout](#))(libct_context_t *context, libct_device_t *device)
Function pointer to the application callback to receive timed out notification in response to calling [libct_connect\(\)](#).
- void(* [on_device_disconnected](#))(libct_context_t *context, libct_device_t *device)
Function pointer to the application callback to receive disconnect notification.
- void(* [on_start_monitoring](#))(libct_context_t *context, libct_device_t *device, int status)
Function pointer to the application callback to receive notification in response to calling [libct_start_monitoring\(\)](#).
- void(* [on_stop_monitoring](#))(libct_context_t *context, libct_device_t *device, int status)
Function pointer to the application callback to receive notification in response to calling [libct_stop_monitoring\(\)](#).
- void(* [on_start_measuring](#))(libct_context_t *context, libct_device_t *device, int status)
Function pointer to the application callback to receive notification in response to calling [libct_start_measuring\(\)](#).
- void(* [on_stop_measuring](#))(libct_context_t *context, libct_device_t *device, int status)
Function pointer to the application callback to receive notification in response to calling [libct_stop_measuring\(\)](#).
- void(* [on_data_received](#))(libct_context_t *context, libct_device_t *device, libct_stream_data_t *data)
Function pointer to the application callback to receive data notifications.
- void(* [on_data_error](#))(libct_context_t *context, libct_device_t *device, const char *error)

- Function pointer to the application callback to receive data error notification.*

 - void(* [on_rd_snr_min_rsp](#))(libct_context_t *context, libct_device_t *device, int snr, int status)
- Function pointer to the application callback to receive status in response to calling [libct_rd_snr_min\(\)](#).*

 - void(* [on_wrt_snr_min_rsp](#))(libct_context_t *context, libct_device_t *device, int status)
- Function pointer to the application callback to receive status in response to calling [libct_wrt_snr_min\(\)](#).*

 - void(* [on_rd_display_state_rsp](#))(libct_context_t *context, libct_device_t *device, unsigned char state, int status)
- Function pointer to the application callback to receive status in response to calling [libct_rd_display_state\(\)](#).*

 - void(* [on_wrt_display_state_rsp](#))(libct_context_t *context, libct_device_t *device, int status)
- Function pointer to the application callback to receive status in response to calling [libct_wrt_display_state\(\)](#).*

 - void(* [on_rd_recal_itvl_rsp](#))(libct_context_t *context, libct_device_t *device, unsigned int itvl, int status)
- Function pointer to the application callback to receive status in response to calling [libct_rd_recal_itvl\(\)](#).*

 - void(* [on_wrt_recal_itvl_rsp](#))(libct_context_t *context, libct_device_t *device, int status)
- Function pointer to the application callback to receive status in response to calling [libct_wrt_recal_itvl\(\)](#).*

 - void(* [on_rd_cuff_pressure_rsp](#))(libct_context_t *context, libct_device_t *device, int pressure, int status)
- Function pointer to the application callback to receive status in response to calling [libct_rd_cuff_pressure\(\)](#).*

 - void(* [on_vent_cuff_rsp](#))(libct_context_t *context, libct_device_t *device, int status)
- Function pointer to the application callback to receive status in response to calling [libct_vent_cuff\(\)](#).*

 - void(* [on_clr_status_rsp](#))(libct_context_t *context, libct_device_t *device, int status)
- Function pointer to the application callback to receive status in response to calling [libct_clr_status\(\)](#).*

 - void(* [on_diag_flush_rsp](#))(libct_context_t *context, libct_device_t *device, int status)
- Function pointer to the application callback to receive status in response to calling [libct_diag_flush\(\)](#).*

 - void(* [on_wrt_waveform_clamping](#))(libct_context_t *context, libct_device_t *device, int status)
- Function pointer to the application callback to receive status in response to calling [libct_wrt_waveform_clamping\(\)](#).*

 - void(* [on_rd_waveform_clamping](#))(libct_context_t *context, libct_device_t *device, unsigned char value, int status)
- Function pointer to the application callback to receive status in response to calling [libct_rd_waveform_clamping\(\)](#).*

 - void(* [on_rd_vitals_filter](#))(libct_context_t *context, libct_device_t *device, int value, int status)
- Function pointer to the application callback to receive status in response to calling [libct_rd_vitals_filter\(\)](#).*

 - void(* [on_wrt_vitals_filter](#))(libct_context_t *context, libct_device_t *device, int status)
- Function pointer to the application callback to receive status in response to calling [libct_wrt_vitals_filter\(\)](#).*

 - void(* [on_rd_motion_timeout](#))(libct_context_t *context, libct_device_t *device, int timeout, int status)
- Function pointer to the application callback to receive status in response to calling [libct_rd_motion_timeout\(\)](#).*

 - void(* [on_rd_persistent_log](#))(libct_context_t *context, libct_device_t *device, const char *log, unsigned int len, int status)
- Function pointer to the application callback to receive status in response to calling [libct_rd_persistent_log\(\)](#).*

7.1.1 Detailed Description

Structure used to provide asynchronous notifications to the application.

This structure is a container of function pointers to your application callback functions to receive asynchronous notifications. Note you are not required to implement all callback functions. Instead, initialize the [libct_app_callbacks_t](#) object to zeros and then set only the function pointers to the callback functions you care about. However, you must implement at least the following callbacks to connect and receive data from the device.

- [on_device_discovered\(\)](#)
- [on_device_connected_ready\(\)](#)
- [on_device_disconnected\(\)](#)

- [on_data_received\(\)](#)

IMPORTANT: When implementing a callback function, you must include **LIBCTAPI** in the function signature to specify the calling convention. This ensures the application and library are using the same calling convention to prevent corrupting the stack. Some platforms, such as Windows, support many calling conventions and **LIBCTAPI** will be set to the default one. If you don't specify **LIBCTAPI** in the callback implementation, the application source code may not compile or serious failures may occur at runtime due to stack corruption. See the sample implementations included with the member descriptions below for details.

7.1.2 Field Documentation

7.1.2.1 on_device_discovered

```
void( * libct_app_callbacks_t::on_device_discovered) (libct_context_t *context, libct_device_t *device)
```

Function pointer to the application callback to receive scan notifications in response to calling [libct_start_discovery\(\)](#).

These notifications are sent to the application during device discovery to notify a discovered device when scanning for Caretaker devices.

The following example illustrates a sample implementation of the callback. Note you must specify **LIBCTAPI** in the function signature.

```
void LIBCTAPI on_device_discovered(libct_context_t* context,
    libct_device_t* device)
{
    const char* address = libct_device_get_address(device);
    const char* name = libct_device_get_name(device);

    // do something
}
```

And you can set the function pointer as follows.

```
libct_app_callbacks_t callbacks = {0};
callbacks.on_device_discovered = on_device_discovered;
```

Parameters

<i>context</i>	The context returned from libct_init() .
<i>device</i>	<p>The discovered device.</p> <p>Note: This device object will be recycled when the callback returns so do not save the device pointer. Instead, make a copy of the device info if needed. Also, since a connection is not established to the device at this point, device functions that require the connection to be established will not return anything useful. You can only call the following functions safely on the device object passed to this callback.</p> <ul style="list-style-type: none">• libct_device_t::get_address()• libct_device_t::get_name()

See also

[on_discovery_timeout\(\)](#)
[on_discovery_failed\(\)](#)

7.1.2.2 on_discovery_timeout

```
void( * libct_app_callbacks_t::on_discovery_timeout) (libct_context_t *context)
```

Function pointer to the application callback to receive timeout notification in response to calling [libct_start_discovery\(\)](#).

The following example illustrates a sample implementation of the callback. Note you must specify **LIBCTAPI** in the function signature.

```
void LIBCTAPI on_discovery_timeout(libct_context_t* context)
{
    // do something
}
```

And you can set the function pointer as follows.

```
libct_app_callbacks_t callbacks = {0};
callbacks.on_discovery_timeout = on_discovery_timeout;
```

Parameters

<i>context</i>	The context returned from libct_init() .
----------------	--

See also

[on_device_discovered\(\)](#)
[on_discovery_failed\(\)](#)

7.1.2.3 on_discovery_failed

```
void( * libct_app_callbacks_t::on_discovery_failed) (libct_context_t *context, int error)
```

Function pointer to the application callback to receive error notification in response to calling [libct_start_discovery\(\)](#).

The following example illustrates a sample implementation of the callback. Note you must specify **LIBCTAPI** in the function signature.

```
void LIBCTAPI on_discovery_failed(libct_context_t* context, int error)
{
    // do something
}
```

And you can set the function pointer as follows.

```
libct_app_callbacks_t callbacks = {0};
callbacks.on_discovery_failed = on_discovery_failed;
```

Parameters

<i>context</i>	The context returned from libct_init() .
<i>error</i>	Generic error code describing the failure.

See also

[on_device_discovered\(\)](#)
[on_discovery_timedout\(\)](#)

7.1.2.4 on_device_connected_not_ready

```
void( * libct_app_callbacks_t::on_device_connected_not_ready) (libct_context_t *context, libct_device_t *device)
```

Function pointer to the application callback to receive early connection notification in response to calling [libct_connect\(\)](#).

This is an early notification to allow the application to update the device connection status, however, the device is not ready for IO at this stage of the connection sequence.

The following example illustrates a sample implementation of the callback. Note you must specify **LIBCTAPI** in the function signature.

```
void LIBCTAPI on_device_connected_not_ready(  
    libct_context_t* context, libct_device_t* device)  
{  
    // update connection status  
}
```

And you can set the function pointer as follows.

```
libct_app_callbacks_t callbacks = {0};  
callbacks.on_device_connected_not_ready =  
    on_device_connected_not_ready;
```

Parameters

<i>context</i>	The context returned from libct_init() .
<i>device</i>	The connected device. Note: Since the connection is established to the device at this point there is no restriction on which device functions you can call to obtain information about the device.

See also

[on_connect_error\(\)](#)
[on_connect_timedout\(\)](#)

7.1.2.5 on_device_connected_ready

```
void( * libct_app_callbacks_t::on_device_connected_ready) (libct_context_t *context, libct_device_t *device)
```

Function pointer to the application callback to receive connection notification in response to calling [libct_connect\(\)](#).

At this stage, the device is ready for IO and the application can issue requests to the device.

The following example illustrates a sample implementation of the callback. Note you must specify **LIBCTAPI** in the function signature.

```
void LIBCTAPI on_device_connected_ready(libct_context_t* context,
    libct_device_t* device)
{
    // do something
}
```

And you can set the function pointer as follows.

```
libct_app_callbacks_t callbacks = {0};
callbacks.on_device_connected_ready =
    on_device_connected_ready;
```

Parameters

<i>context</i>	The context returned from libct_init() .
<i>device</i>	The connected device. Note: Since the connection is established to the device at this point there is no restriction on which device functions you can call to obtain information about the device.

See also

[on_connect_error\(\)](#)
[on_connect_timedout\(\)](#)

7.1.2.6 on_connect_error

```
void( * libct_app_callbacks_t::on_connect_error) (libct_context_t *context, libct_device_t *device, const char *error)
```

Function pointer to the application callback to receive error notification in response to calling [libct_connect\(\)](#).

The following example illustrates a sample implementation of the callback. Note you must specify **LIBCTAPI** in the function signature..

```
void LIBCTAPI on_connect_error(libct_context_t* context,
    libct_device_t* device, const char* error)
{
    // do something
}
```

And you can set the function pointer as follows.

```
libct_app_callbacks_t callbacks = {0};
callbacks.on_connect_error = on_connect_error;
```

Parameters

<i>context</i>	The context returned from libct_init() .
<i>device</i>	The affected device. Note: Since there is no connection to the device the device functions requiring a connection will not return anything useful.
<i>error</i>	String describing error.

See also

[on_device_connected_not_ready\(\)](#)
[on_device_connected_ready\(\)](#)
[on_connect_timedout\(\)](#)

7.1.2.7 on_connect_timedout

```
void( * libct_app_callbacks_t::on_connect_timedout) (libct_context_t *context, libct_device_t *device)
```

Function pointer to the application callback to receive timed out notification in response to calling [libct_connect\(\)](#).

The following example illustrates a sample implementation of the callback. Note you must specify **LIBCTAPI** in the function signature.

```
void LIBCTAPI on_connect_timedout(libct_context_t* context,  
    libct_device_t* device)  
{  
    // do something  
}
```

And you can set the function pointer as follows.

```
libct_app_callbacks_t callbacks = {0};  
callbacks.on_connect_timedout = on_connect_timedout;
```

Parameters

<i>context</i>	The context returned from libct_init() .
<i>device</i>	The affected device. Note: Since there is no connection to the device the device functions requiring a connection will not return anything useful.

See also

[on_device_connected_not_ready\(\)](#)
[on_device_connected_ready\(\)](#)
[on_connect_error\(\)](#)

7.1.2.8 on_device_disconnected

```
void( * libct_app_callbacks_t::on_device_disconnected) (libct_context_t *context, libct_device_t *device)
```

Function pointer to the application callback to receive disconnect notification.

The disconnect notification is notified after the connection is established with the device and the connection is lost such as when the device moves out of range and disconnects.

The disconnect notification is also notified when the application calls `libct_disconnect()` to disconnect explicitly. However, the notification to the application is not guaranteed to occur for this scenario, which should be okay since the application initiated the disconnect.

IMPORTANT: The application must not call `libct_disconnect()` or `libct_deinit()` from within this or any library callback function. Callbacks are called from internal library threads that these functions attempt to kill. As such, `libct_disconnect()` and `libct_deinit()` must only be called from an application thread.

The following example illustrates a sample implementation of the callback. Note you must specify **LIBCTAPI** in the function signature.

```
void LIBCTAPI on_device_disconnected(libct_context_t* context,
    libct_device_t* device)
{
    // do something
}
```

And you can set the function pointer as follows.

```
libct_app_callbacks_t callbacks = {0};
callbacks.on_device_disconnected = on_device_disconnected;
```

Parameters

<i>context</i>	The context returned from <code>libct_init()</code> .
<i>device</i>	The disconnected device.

7.1.2.9 on_start_monitoring

```
void( * libct_app_callbacks_t::on_start_monitoring) (libct_context_t *context, libct_device_t *device, int status)
```

Function pointer to the application callback to receive notification in response to calling `libct_start_monitoring()`.

This notification will be sent only once (one-shot) to notify success or failure after calling `libct_start_monitoring()`. If success, the application `on_data_received()` will be notified repeatedly with data from the device until `libct_stop_monitoring()` is called subsequently or the device is disconnected.

The following example illustrates a sample implementation of the callback. Note you must specify **LIBCTAPI** in the function signature.

```
void LIBCTAPI on_start_monitoring(libct_context_t* context,
    libct_device_t* device, int status)
{
    // do something
}
```

And you can set the function pointer as follows.

```
libct_app_callbacks_t callbacks = {0};
callbacks.on_start_monitoring = on_start_monitoring;
```

Parameters

<i>context</i>	The context returned from libct_init() .
<i>device</i>	The device being monitored.
<i>status</i>	Status indicating success or failure: zero on success and non-zero otherwise.

See also

[on_data_received\(\)](#)

7.1.2.10 on_stop_monitoring

```
void( * libct_app_callbacks_t::on_stop_monitoring) (libct_context_t *context, libct_device_t *device, int status)
```

Function pointer to the application callback to receive notification in response to calling [libct_stop_monitoring\(\)](#).

This notification will be sent only once (one-shot) to notify success or failure after calling [libct_stop_monitoring\(\)](#). If success, the application [on_data_received\(\)](#) will stop receiving data notifications.

The following example illustrates a sample implementation of the callback. Note you must specify **LIBCTAPI** in the function signature.

```
void LIBCTAPI on_stop_monitoring(libct_context_t* context,
    libct_device_t* device, int status)
{
    // do something
}
```

And you can set the function pointer as follows.

```
libct_app_callbacks_t callbacks = {0};
callbacks.on_stop_monitoring = on_stop_monitoring;
```

Parameters

<i>context</i>	The context returned from libct_init() .
<i>device</i>	The device being monitored.
<i>status</i>	Status indicating success or failure: zero on success and non-zero otherwise.

7.1.2.11 on_start_measuring

```
void( * libct_app_callbacks_t::on_start_measuring) (libct_context_t *context, libct_device_t *device, int status)
```

Function pointer to the application callback to receive notification in response to calling [libct_start_measuring\(\)](#).

This notification will be sent only once (one-shot) to notify success or failure after calling [libct_start_measuring\(\)](#). If success, the device will begin pulse decomposition analysis (PDA) and vital sign measurements (blood pressure, heart rate, etc) will be notified to the application [on_data_received\(\)](#) callback.

The following example illustrates a sample implementation of the callback. Note you must specify **LIBCTAPI** in the function signature.

```
void LIBCTAPI on_start_measuring(libct_context_t* context,
                                libct_device_t* device, int status)
{
    // do something
}
```

And you can set the function pointer as follows.

```
libct_app_callbacks_t callbacks = {0};
callbacks.on_start_measuring = on_start_measuring;
```

Parameters

<i>context</i>	The context returned from libct_init() .
<i>device</i>	The device providing measurements.
<i>status</i>	Status indicating success or failure: zero on success and non-zero otherwise.

7.1.2.12 on_stop_measuring

```
void( * libct_app_callbacks_t::on_stop_measuring) (libct_context_t *context, libct_device_t *device, int status)
```

Function pointer to the application callback to receive notification in response to calling [libct_stop_measuring\(\)](#).

This notification will be sent only once (one-shot) to notify success or failure after calling [libct_stop_measuring\(\)](#). If success, the device will stop pulse decomposition analysis (PDA) and vital sign measurements (blood pressure, heart rate, etc) will stop.

The following example illustrates a sample implementation of the callback. Note you must specify **LIBCTAPI** in the function signature.

```
void LIBCTAPI on_stop_measuring(libct_context_t* context,
                                libct_device_t* device, int status)
{
    // do something
}
```

And you can set the function pointer as follows.

```
libct_app_callbacks_t callbacks = {0};
callbacks.on_stop_measuring = on_stop_measuring;
```

Parameters

<i>context</i>	The context returned from libct_init() .
<i>device</i>	The device providing measurements.
<i>status</i>	Status indicating success or failure: zero on success and non-zero otherwise.

7.1.2.13 on_data_received

```
void( * libct_app_callbacks_t::on_data_received) (libct_context_t *context, libct_device_t *device, libct_stream_data_t *data)
```

Function pointer to the application callback to receive data notifications.

These notifications are sent repeatedly to the application to hand-off data received from the device some time after calling [libct_start_monitoring\(\)](#) successfully.

Data notified via this callback depends on the [monitor flags](#) passed to [libct_start_monitoring\(\)](#) and whether or not [libct_start_measuring\(\)](#) was called to start taking vital sign measurements. See [libct_stream_data_t](#) for data details.

The following example illustrates a sample implementation of the callback. Note you must specify **LIBCTAPI** in the function signature.

```
void LIBCTAPI on_data_received(libct_context_t* context,
    libct_device_t* device, libct_stream_data_t* data)
{
    // handle data
}
```

And you can set the function pointer as follows.

```
libct_app_callbacks_t callbacks = {0};
callbacks.on_data_received = on_data_received;
```

Parameters

<i>context</i>	The context returned from libct_init() .
<i>device</i>	The device originating the data.
<i>data</i>	Stream packet containing the data received from the device. NOTE: The stream data packet is created with dynamic memory that will be freed after the callback returns. So you should not save pointer(s) to the data, instead copy individual fields into application memory as needed if you need to access it after on_data_received() returns. Do not copy the entire libct_stream_data_t structure as it is a structure of pointers and doing so will be saving pointers to freed memory after the callback returned.

See also

[on_data_error\(\)](#)

7.1.2.14 on_data_error

```
void( * libct_app_callbacks_t::on_data_error) (libct_context_t *context, libct_device_t *device,
const char *error)
```

Function pointer to the application callback to receive data error notification.

This notification is sent if the library encounters error receiving or processing data.

The following example illustrates a sample implementation of the callback. Note you must specify **LIBCTAPI** in the function signature.

```
void LIBCTAPI on_data_error(libct_context_t* context,
    libct_device_t* device, const char* error)
{
    // handle error
}
```

And you can set the function pointer as follows.

```
libct_app_callbacks_t callbacks = {0};
callbacks.on_data_error = on_data_error;
```

Parameters

<i>context</i>	The context returned from libct_init() .
<i>device</i>	The affected device.
<i>error</i>	String describing the error.

7.1.2.15 on_rd_snr_min_rsp

```
void( * libct_app_callbacks_t::on_rd_snr_min_rsp) (libct_context_t *context, libct_device_t *device, int snr, int status)
```

Function pointer to the application callback to receive status in response to calling [libct_rd_snr_min\(\)](#).

This notification will be sent only once (one-shot) to notify success or failure after calling [libct_rd_snr_min\(\)](#).

The following example illustrates a sample implementation of the callback. Note you must specify **LIBCTAPI** in the function signature.

```
void LIBCTAPI on_rd_snr_min_rsp(libct_context_t* context,
    libct_device_t* device, int snr, int status)
{
    // do something
}
```

And you can set the function pointer as follows.

```
libct_app_callbacks_t callbacks = {0};  
callbacks.on_rd_snr_min_rsp = on_rd_snr_min_rsp;
```

Parameters

<i>context</i>	The context returned from libct_init() .
<i>device</i>	The device associated with the context.
<i>snr</i>	Minimum signal-to-noise value on success.
<i>status</i>	Status indicating success or failure: zero on success and non-zero otherwise.

7.1.2.16 on_wrt_snr_min_rsp

```
void( * libct_app_callbacks_t::on_wrt_snr_min_rsp) (libct_context_t *context, libct_device_t *device, int status)
```

Function pointer to the application callback to receive status in response to calling [libct_wrt_snr_min\(\)](#).

This notification will be sent only once (one-shot) to notify success or failure after calling [libct_wrt_snr_min\(\)](#).

The following example illustrates a sample implementation of the callback. Note you must specify **LIBCTAPI** in the function signature.

```
void LIBCTAPI on_wrt_snr_min_rsp(libct_context_t* context,  
    libct_device_t* device, int status)  
{  
    // do something  
}
```

And you can set the function pointer as follows.

```
libct_app_callbacks_t callbacks = {0};  
callbacks.on_wrt_snr_min_rsp = on_wrt_snr_min_rsp;
```

Parameters

<i>context</i>	The context returned from libct_init() .
<i>device</i>	The device associated with the context.
<i>status</i>	Status indicating success or failure: zero on success and non-zero otherwise.

7.1.2.17 on_rd_display_state_rsp

```
void( * libct_app_callbacks_t::on_rd_display_state_rsp) (libct_context_t *context, libct_device_t *device, unsigned char state, int status)
```

Function pointer to the application callback to receive status in response to calling [libct_rd_display_state\(\)](#).

This notification will be sent only once (one-shot) to notify success or failure after calling [libct_rd_display_state\(\)](#).

The following example illustrates a sample implementation of the callback. Note you must specify **LIBCTAPI** in the function signature.

```
void LIBCTAPI on_rd_display_state_rsp(libct_context_t* context,
    libct_device_t* device, unsigned char state, int status)
{
    // do something
}
```

And you can set the function pointer as follows.

```
libct_app_callbacks_t callbacks = {0};
callbacks.on_rd_display_state_rsp =
    on_rd_display_state_rsp;
```

Parameters

<i>context</i>	The context returned from libct_init() .
<i>device</i>	The device associated with the context.
<i>state</i>	The device display state on success: 0 = off, 1 = on.
<i>status</i>	Status indicating success or failure: zero on success and non-zero otherwise.

7.1.2.18 on_wrt_display_state_rsp

```
void( * libct_app_callbacks_t::on_wrt_display_state_rsp) (libct_context_t *context, libct_↵
device_t *device, int status)
```

Function pointer to the application callback to receive status in response to calling [libct_wrt_display_state\(\)](#).

This notification will be sent only once (one-shot) to notify success or failure after calling [libct_wrt_display_state\(\)](#).

The following example illustrates a sample implementation of the callback. Note you must specify **LIBCTAPI** in the function signature.

```
void LIBCTAPI on_wrt_display_state_rsp(libct_context_t* context,
    libct_device_t* device, int status)
{
    // do something
}
```

And you can set the function pointer as follows.

```
libct_app_callbacks_t callbacks = {0};
callbacks.on_wrt_display_state_rsp =
    on_wrt_display_state_rsp;
```

Parameters

<i>context</i>	The context returned from libct_init() .
<i>device</i>	The device associated with the context.
<i>status</i>	Status indicating success or failure: zero on success and non-zero otherwise.

7.1.2.19 `on_rd_recal_itvl_rsp`

```
void( * libct_app_callbacks_t::on_rd_recal_itvl_rsp) (libct_context_t *context, libct_device_t *device, unsigned int itvl, int status)
```

Function pointer to the application callback to receive status in response to calling [libct_rd_recal_itvl\(\)](#).

This notification will be sent only once (one-shot) to notify success or failure after calling [libct_rd_recal_itvl\(\)](#).

The following example illustrates a sample implementation of the callback. Note you must specify **LIBCTAPI** in the function signature.

```
void LIBCTAPI on_rd_recal_itvl_rsp(libct_context_t* context,
    libct_device_t* device, unsigned int itvl, int status)
{
    // do something
}
```

And you can set the function pointer as follows.

```
libct_app_callbacks_t callbacks = {0};
callbacks.on_rd_recal_itvl_rsp = on_rd_recal_itvl_rsp;
```

Parameters

<i>context</i>	The context returned from libct_init() .
<i>device</i>	The device associated with the context.
<i>itvl</i>	The recalibration interval in minutes on success.
<i>status</i>	Status indicating success or failure: zero on success and non-zero otherwise.

7.1.2.20 `on_wrt_recal_itvl_rsp`

```
void( * libct_app_callbacks_t::on_wrt_recal_itvl_rsp) (libct_context_t *context, libct_device_t *device, int status)
```

Function pointer to the application callback to receive status in response to calling [libct_wrt_recal_itvl\(\)](#).

This notification will be sent only once (one-shot) to notify success or failure after calling [libct_wrt_recal_itvl\(\)](#).

The following example illustrates a sample implementation of the callback. Note you must specify **LIBCTAPI** in the function signature.

```
void LIBCTAPI on_wrt_recal_itvl_rsp(libct_context_t* context,
    libct_device_t* device, int status)
{
    // do something
}
```

And you can set the function pointer as follows.

```
libct_app_callbacks_t callbacks = {0};
callbacks.on_wrt_recal_itvl_rsp = on_wrt_recal_itvl_rsp;
```

Parameters

<i>context</i>	The context returned from libct_init() .
<i>device</i>	The device associated with the context.
<i>status</i>	Status indicating success or failure: zero on success and non-zero otherwise.

7.1.2.21 on_rd_cuff_pressure_rsp

```
void( * libct_app_callbacks_t::on_rd_cuff_pressure_rsp) (libct_context_t *context, libct_↵
device_t *device, int pressure, int status)
```

Function pointer to the application callback to receive status in response to calling [libct_rd_cuff_pressure\(\)](#).

This notification will be sent only once (one-shot) to notify success or failure after calling [libct_rd_cuff_pressure\(\)](#).

The following example illustrates a sample implementation of the callback. Note you must specify **LIBCTAPI** in the function signature.

```
void LIBCTAPI on_rd_cuff_pressure_rsp(libct_context_t* context,
    libct_device_t* device, int status)
{
    // do something
}
```

And you can set the function pointer as follows.

```
libct_app_callbacks_t callbacks = {0};
callbacks.on_rd_cuff_pressure_rsp =
    on_rd_cuff_pressure_rsp;
```

Parameters

<i>context</i>	The context returned from libct_init() .
<i>device</i>	The device associated with the context.
<i>pressure</i>	The cuff pressure in mmHg on success.
<i>status</i>	Status indicating success or failure: zero on success and non-zero otherwise.

7.1.2.22 on_vent_cuff_rsp

```
void( * libct_app_callbacks_t::on_vent_cuff_rsp) (libct_context_t *context, libct_device_t *device, int status)
```

Function pointer to the application callback to receive status in response to calling [libct_vent_cuff\(\)](#).

This notification will be sent only once (one-shot) to notify success or failure after calling [libct_vent_cuff\(\)](#).

The following example illustrates a sample implementation of the callback. Note you must specify **LIBCTAPI** in the function signature.

```
void LIBCTAPI on_vent_cuff_rsp(libct_context_t* context,
    libct_device_t* device, int status)
{
    // do something
}
```

And you can set the function pointer as follows.

```
libct_app_callbacks_t callbacks = {0};
callbacks.on_vent_cuff_rsp = on_vent_cuff_rsp;
```

Parameters

<i>context</i>	The context returned from libct_init() .
<i>device</i>	The device associated with the context.
<i>status</i>	Status indicating success or failure: zero on success and non-zero otherwise.

7.1.2.23 on_clr_status_rsp

```
void( * libct_app_callbacks_t::on_clr_status_rsp) (libct_context_t *context, libct_device_t *device, int status)
```

Function pointer to the application callback to receive status in response to calling [libct_clr_status\(\)](#).

This notification will be sent only once (one-shot) to notify success or failure after calling [libct_clr_status\(\)](#).

The following example illustrates a sample implementation of the callback. Note you must specify **LIBCTAPI** in the function signature.

```
void LIBCTAPI on_clr_status_rsp(libct_context_t* context,
    libct_device_t* device, int status)
{
    // do something
}
```

And you can set the function pointer as follows.

```
libct_app_callbacks_t callbacks = {0};
callbacks.on_clr_status_rsp = on_clr_status_rsp;
```


Parameters

<i>context</i>	The context returned from libct_init() .
<i>device</i>	The device associated with the context.
<i>status</i>	Status indicating success or failure: zero on success and non-zero otherwise.

7.1.2.24 `on_diag_flush_rsp`

```
void( * libct_app_callbacks_t::on_diag_flush_rsp) (libct_context_t *context, libct_device_t *device, int status)
```

Function pointer to the application callback to receive status in response to calling [libct_diag_flush\(\)](#).

This notification will be sent only once (one-shot) to notify success or failure after calling [libct_diag_flush\(\)](#).

The following example illustrates a sample implementation of the callback. Note you must specify **LIBCTAPI** in the function signature.

```
void LIBCTAPI on_diag_flush_rsp(libct_context_t* context,
    libct_device_t* device, int status)
{
    // do something
}
```

And you can set the function pointer as follows.

```
libct_app_callbacks_t callbacks = {0};
callbacks.on_diag_flush_rsp = on_diag_flush_rsp;
```

Parameters

<i>context</i>	The context returned from libct_init() .
<i>device</i>	The device associated with the context.
<i>status</i>	Status indicating success or failure: zero on success and non-zero otherwise.

7.1.2.25 `on_wrt_waveform_clamping`

```
void( * libct_app_callbacks_t::on_wrt_waveform_clamping) (libct_context_t *context, libct_device_t *device, int status)
```

Function pointer to the application callback to receive status in response to calling [libct_wrt_waveform_clamping\(\)](#).

This notification will be sent only once (one-shot) to notify success or failure after calling [libct_wrt_waveform_clamping\(\)](#).

The following example illustrates a sample implementation of the callback. Note you must specify **LIBCTAPI** in the function signature.

```
void LIBCTAPI on_wrt_waveform_clamping(libct_context_t* context,
                                      libct_device_t* device
                                      int status)
{
    // do something
}
```

And you can set the function pointer as follows.

```
libct_app_callbacks_t callbacks = {0};
callbacks.on_wrt_waveform_clamping =
    on_wrt_waveform_clamping;
```

Parameters

<i>context</i>	The context returned from libct_init() .
<i>device</i>	The device associated with the context.
<i>status</i>	Status indicating success or failure: zero on success and non-zero otherwise.

7.1.2.26 on_rd_waveform_clamping

```
void( * libct_app_callbacks_t::on_rd_waveform_clamping) (libct_context_t *context, libct_↵
device_t *device, unsigned char value, int status)
```

Function pointer to the application callback to receive status in response to calling [libct_rd_waveform_clamping\(\)](#).

This notification will be sent only once (one-shot) to notify success or failure after calling [libct_wrt_waveform_↵
clamping\(\)](#).

The following example illustrates a sample implementation of the callback. Note you must specify **LIBCTAPI** in the function signature.

```
void LIBCTAPI on_rd_waveform_clamping(libct_context_t* context,
                                      libct_device_t* device,
                                      unsigned char value
                                      int status)
{
    // do something
}
```

And you can set the function pointer as follows.

```
libct_app_callbacks_t callbacks = {0};
callbacks.on_rd_waveform_clamping =
    on_rd_waveform_clamping;
```

Parameters

<i>context</i>	The context returned from libct_init() .
<i>device</i>	The device associated with the context.
<i>value</i>	Clamp setting: 1 = ON, 0 = OFF
<i>status</i>	Status indicating success or failure: zero on success and non-zero otherwise.

7.1.2.27 on_rd_vitals_filter

```
void( * libct_app_callbacks_t::on_rd_vitals_filter) (libct_context_t *context, libct_device_t *device, int value, int status)
```

Function pointer to the application callback to receive status in response to calling [libct_rd_vitals_filter\(\)](#).

This notification will be sent only once (one-shot) to notify success or failure after calling [libct_rd_vitals_filter\(\)](#).

The following example illustrates a sample implementation of the callback. Note you must specify **LIBCTAPI** in the function signature.

```
void LIBCTAPI on_rd_vitals_filter(libct_context_t* context,
                                libct_device_t* device, int status)
{
    // do something
}
```

And you can set the function pointer as follows.

```
libct_app_callbacks_t callbacks = {0};
callbacks.on_rd_vitals_filter = on_rd_vitals_filter;
```

Parameters

<i>context</i>	The context returned from libct_init() .
<i>device</i>	The device associated with the context.
<i>value</i>	The median filter value when success. 0 = Disabled, 1 = Enabled.
<i>status</i>	Status indicating success or failure: zero on success and nonzero otherwise.

7.1.2.28 on_wrt_vitals_filter

```
void( * libct_app_callbacks_t::on_wrt_vitals_filter) (libct_context_t *context, libct_device_t *device, int status)
```

Function pointer to the application callback to receive status in response to calling [libct_wrt_vitals_filter\(\)](#).

This notification will be sent only once (one-shot) to notify success or failure after calling [libct_wrt_vitals_filter\(\)](#).

The following example illustrates a sample implementation of the callback. Note you must specify **LIBCTAPI** in the function signature.

```
void LIBCTAPI on_wrt_vitals_filter(libct_context_t* context,
                                libct_device_t* device, int status)
{
    // do something
}
```

And you can set the function pointer as follows.

```
libct_app_callbacks_t callbacks = {0};
callbacks.on_wrt_vitals_filter = on_wrt_vitals_filter;
```

Parameters

<i>context</i>	The context returned from libct_init() .
<i>device</i>	The device associated with the context.
<i>status</i>	Status indicating success or failure: zero on success and non-zero otherwise.

7.1.2.29 on_rd_motion_timeout

```
void( * libct_app_callbacks_t::on_rd_motion_timeout) (libct_context_t *context, libct_device_t *device, int timeout, int status)
```

Function pointer to the application callback to receive status in response to calling [libct_rd_motion_timeout\(\)](#).

This notification will be sent only once (one-shot) to notify success or failure after calling [libct_rd_motion_timeout\(\)](#).

The following example illustrates a sample implementation of the callback. Note you must specify **LIBCTAPI** in the function signature.

```
void LIBCTAPI on_rd_motion_timeout(libct_context_t* context,
    libct_device_t* device, int status, int timeout)
{
    // do something
}
```

And you can set the function pointer as follows.

```
libct_app_callbacks_t callbacks = {0};
callbacks.on_rd_motion_timeout = on_rd_motion_timeout;
```

Parameters

<i>context</i>	The context returned from libct_init() .
<i>device</i>	The device associated with the context.
<i>status</i>	Status indicating success or failure: zero on success and non-zero otherwise.
<i>timeout</i>	Motion timeout value in seconds.

7.1.2.30 on_rd_persistent_log

```
void( * libct_app_callbacks_t::on_rd_persistent_log) (libct_context_t *context, libct_device_t *device, const char *log, unsigned int len, int status)
```

Function pointer to the application callback to receive status in response to calling [libct_rd_persistent_log\(\)](#).

This notification will be sent only once (one-shot) to notify success or failure after calling [libct_rd_persistent_log\(\)](#).

The following example illustrates a sample implementation of the callback. Note you must specify **LIBCTAPI** in the function signature.

```
void LIBCTAPI on_rd_persistent_log(libct_context_t* context,
                                  libct_device_t* device,
                                  const char* log,
                                  unsigned int len,
                                  int timeout)
{
    // do something
}
```

And you can set the function pointer as follows.

```
libct_app_callbacks_t callbacks = {0};
callbacks.on_rd_persistent_log = on_rd_persistent_log;
```

Parameters

<i>context</i>	The context returned from libct_init() .
<i>device</i>	The device associated with the context.
<i>log</i>	The device persistent log.
<i>len</i>	Log length.
<i>status</i>	Status indicating success or failure: zero on success and non-zero otherwise.

The documentation for this struct was generated from the following file:

- caretaker.h

7.2 libct_battery_info_t Class Reference

Battery info data point within the [libct_stream_data_t](#) packet.

Data Fields

- bool [valid](#)
The other fields are valid when this field is non-zero (true) and invalid otherwise.
- int [voltage](#)
The battery voltage in millivolts.
- unsigned int [timestamp](#)
Time stamp from the device associated with the data.

7.2.1 Detailed Description

Battery info data point within the [libct_stream_data_t](#) packet.

7.2.2 Field Documentation

7.2.2.1 voltage

```
int libct_battery_info_t::voltage
```

The battery voltage in millivolts.

7.2.2.2 timestamp

```
unsigned int libct_battery_info_t::timestamp
```

Time stamp from the device associated with the data.

The documentation for this class was generated from the following file:

- caretaker.h

7.3 libct_bp_settings_t Class Reference

Structure to write the caretaker manual blood pressure settings.

Data Fields

- unsigned short [systolic](#)
Systolic pressure setting used for blood pressure calibration.
- unsigned short [diastolic](#)
Diastolic pressure setting used for blood pressure calibration.

7.3.1 Detailed Description

Structure to write the caretaker manual blood pressure settings.

7.3.2 Field Documentation

7.3.2.1 systolic

```
unsigned short libct_bp_settings_t::systolic
```

Systolic pressure setting used for blood pressure calibration.

Acceptable range [80, 250].

7.3.2.2 diastolic

```
unsigned short libct_bp_settings_t::diastolic
```

Diastolic pressure setting used for blood pressure calibration.

Acceptable range [50, 150].

The documentation for this class was generated from the following file:

- caretaker.h

7.4 libct_cal_curve_t Class Reference

Calibration curve data point within the [libct_stream_data_t](#) packet.

Data Fields

- bool [valid](#)
The other fields are valid when this field is non-zero (true) and invalid otherwise.
- int [data_id](#)
Data ID.
- float [val1](#)
Value 1.
- float [val2](#)
Value 2.
- float [val3](#)
Value 3.
- char * [alternateData](#)
Alternate Data - does not come from device.

7.4.1 Detailed Description

Calibration curve data point within the [libct_stream_data_t](#) packet.

Note

The cal curve data is for internal use or research only

7.4.2 Field Documentation

7.4.2.1 valid

```
bool libct_cal_curve_t::valid
```

The other fields are valid when this field is non-zero (true) and invalid otherwise.

7.4.2.2 data_id

```
int libct_cal_curve_t::data_id
```

Data ID.

7.4.2.3 val1

```
float libct_cal_curve_t::val1
```

Value 1.

7.4.2.4 val2

```
float libct_cal_curve_t::val2
```

Value 2.

7.4.2.5 val3

```
float libct_cal_curve_t::val3
```

Value 3.

The documentation for this class was generated from the following file:

- caretaker.h

7.5 libct_cal_t Struct Reference

Structure used to pass calibration data to [libct_start_measuring\(\)](#).

Data Fields

- int [type](#)
Calibration type.
- union {
 - struct {
 - short [posture](#)
Patient posture.
 - } [auto_cal](#)
Calibration configuration when [type](#) is LIBCT_AUTO_CAL.
 - struct {
 - [libct_bp_settings_t](#) **settings**
 - } [manual_cal](#)
Calibration configuration when [type](#) is LIBCT_MANUAL_CAL.
- } [config](#)

Calibration data.

7.5.1 Detailed Description

Structure used to pass calibration data to [libct_start_measuring\(\)](#).

7.5.2 Field Documentation

7.5.2.1 type

```
int libct_cal_t::type
```

Calibration type.

Set to one of the [calibration types](#).

7.5.2.2 posture

```
short libct_cal_t::posture
```

Patient posture.

Set to one of the [patient postures](#).

7.5.2.3 config

```
union { ... } libct_cal_t::config
```

Calibration data.

The documentation for this struct was generated from the following file:

- caretaker.h

7.6 libct_cal_type_t Class Reference

The Caretaker calibration types.

7.6.1 Detailed Description

The Caretaker calibration types.

The documentation for this class was generated from the following file:

- caretaker.h

7.7 libct_context_t Class Reference

An opaque type representing a library instance associated with (or bound to) a device the application is monitoring.

7.7.1 Detailed Description

An opaque type representing a library instance associated with (or bound to) a device the application is monitoring.

The context is used internally to manage the library instance so its data structure is not exposed to the application. As such, the application cannot create a library context explicitly. A library context can only be created by calling [libct_init\(\)](#) to initialize a library instance, which sets the context pointer passed in the first argument. If the call succeeded, the application can use the context to call other library functions, but must call [libct_deinit\(\)](#) to destroy the context when it is no longer needed. Destroying the context releases resources that were allocated when the context was initialized, so the application is required to call [libct_deinit\(\)](#) to release the context, and not doing so will leak system resources.

```
// Initialize library instance, which returns a device context pointer.
libct_context_t* context = NULL;
int status = libct_init(&context, &init_data, &app_callbacks);
if ( LIBCT_FAILED(status) ) {
    // Handle error
    return status;
}

// Connect to a device and monitor data (code not shown)

// Destroy context
libct_deinit(context);
```

The documentation for this class was generated from the following file:

- caretaker.h

7.8 libct_cuff_pressure_t Class Reference

Cuff pressure data point within the [libct_stream_data_t](#) packet.

Data Fields

- bool [valid](#)
The other fields are valid when this field is non-zero (true) and invalid otherwise.
- int [value](#)
cuff pressure actual value.
- int [target](#)
cuff pressure target value.
- int [snr](#)
signal to noise ratio.
- unsigned int [timestamp](#)
Time stamp from the device associated with the data.

7.8.1 Detailed Description

Cuff pressure data point within the [libct_stream_data_t](#) packet.

7.8.2 Field Documentation

7.8.2.1 [valid](#)

```
bool libct_cuff_pressure_t::valid
```

The other fields are valid when this field is non-zero (true) and invalid otherwise.

7.8.2.2 [value](#)

```
int libct_cuff_pressure_t::value
```

cuff pressure actual value.

7.8.2.3 [target](#)

```
int libct_cuff_pressure_t::target
```

cuff pressure target value.

7.8.2.4 snr

```
int libct_cuff_pressure_t::snr
```

signal to noise ratio.

7.8.2.5 timestamp

```
unsigned int libct_cuff_pressure_t::timestamp
```

Time stamp from the device associated with the data.

The documentation for this class was generated from the following file:

- caretaker.h

7.9 libct_device_class_t Class Reference

Classes of devices that can be monitored by this library.

7.9.1 Detailed Description

Classes of devices that can be monitored by this library.

The documentation for this class was generated from the following file:

- caretaker.h

7.10 libct_device_state_t Class Reference

The Caretaker device states.

7.10.1 Detailed Description

The Caretaker device states.

The documentation for this class was generated from the following file:

- caretaker.h

7.11 libct_device_status_t Class Reference

Device status data point within the [libct_stream_data_t](#) packet.

Data Fields

- bool [valid](#)
The other fields are valid when this field is non-zero (true) and invalid otherwise.
- long long [value](#)
Integer value representing logically OR of all status flags, which essentially is the raw value from the device.
- bool [pda_enabled](#)
An indicator of whether the system PDA measurement system is enabled.
- bool [simulation_enabled](#)
An indicator of whether the system is in simulation mode.
- bool [pressure_control_indicator](#)
An indicator of whether the system is currently running closed loop pressure control.
- bool [inflated_indicator](#)
An indicator of whether the system has been inflated to pressure.
- bool [clock_wrap_around](#)
The system clock (time since reset) has wrapped around its index.
- bool [battery_voltage_low](#)
The battery voltage sensor has indicated the battery is near drop-out.
- bool [critical_temperature](#)
The on-board temperature sensor has detected critically high temperature.
- bool [pump_overrun](#)
The pump has violated an overrun condition.
- bool [ble_temperature_sensor_paired](#)
The BLE Temperature Sensor is paired and actively communicating with CareTaker.
- bool [ble_handheld_paired](#)
The a BLE hand-held device is paired and actively communicating with CareTaker.
- bool [ble_stream_control](#)
The current stream control status bit of the BLE stream.
- bool [cellular_control](#)
The current stream control status bit of the cellular stream.
- bool [serial_stream_control](#)
The current stream control status bit of the serial stream.
- bool [auto_cal_mode](#)
The system has been started and running in auto-calibration mode.
- bool [manual_cal_mode](#)
The system has been started and running in manual calibration mode.
- bool [motion_event](#)
The system is having trouble getting a good reading due to too much motion.
- bool [poor_signal](#)
The system failed to calibrate or timed out process signals so measurements were aborted.
- bool [data_valid](#)
There are valid vital signs measurements.
- bool [calibrating](#)
The system is currently calibrating the blood pressure system.
- bool [calibrated](#)

- The system has current valid calibration.*

 - bool [beta_processing](#)

The system has finished finding the oscillometric curve and is processing the beta (offset) value.
 - bool [inflate_failed](#)

Cuff did not inflate to expected value within timeout.
 - bool [calibration_failed](#)

The calibration values were out of range or oscillometric curve had invalid shape.
 - bool [calibration_offset_failed](#)

Too much movement.
 - bool [no_pulse_timeout](#)

The systems has gone greater than 3 minutes without a valid heart beat.
 - bool [cuff_too_loose](#)

The calibration pump up identified the cuff was too loose.
 - bool [cuff_too_tight](#)

The calibration pump up identified the cuff was too tight.
 - bool [weak_signal](#)

Calibration oscillometric curve amplitude is too weak to verify reading.
 - bool [bad_cuff](#)

The cuff is not holding pressure as expected.
 - bool [ble_adv](#)

The bluetooth module is advertising.
 - bool [recal_soon](#)

An automatic recalibration will be occurring shortly.
 - bool [too_many_fails](#)

Auto-calibration failed too many consecutive times try manual calibration.
 - short [autocal_pct](#)

Auto-calibration percentage complete.
 - bool [charging](#)

The device is charging.
 - bool [chargeComplete](#)

Charging complete.
 - short [posture](#)

Posture.

7.11.1 Detailed Description

Device status data point within the [libct_stream_data_t](#) packet.

7.11.2 Field Documentation

7.11.2.1 [simulation_enabled](#)

```
bool libct_device_status_t::simulation_enabled
```

An indicator of whether the system is in simulation mode.

7.11.2.2 `inflated_indicator`

```
bool libct_device_status_t::inflated_indicator
```

An indicator of whether the system has been inflated to pressure.

7.11.2.3 `clock_wrap_around`

```
bool libct_device_status_t::clock_wrap_around
```

The system clock (time since reset) has wrapped around its index.

7.11.2.4 `battery_voltage_low`

```
bool libct_device_status_t::battery_voltage_low
```

The battery voltage sensor has indicated the battery is near drop-out.

7.11.2.5 `pump_overrun`

```
bool libct_device_status_t::pump_overrun
```

The pump has violated an overrun condition.

7.11.2.6 `ble_temperature_sensor_paired`

```
bool libct_device_status_t::ble_temperature_sensor_paired
```

The BLE Temperature Sensor is paired and actively communicating with CareTaker.

7.11.2.7 `ble_stream_control`

```
bool libct_device_status_t::ble_stream_control
```

The current stream control status bit of the BLE stream.

7.11.2.8 manual_cal_mode

```
bool libct_device_status_t::manual_cal_mode
```

The system has been started and running in manual calibration mode.

7.11.2.9 motion_event

```
bool libct_device_status_t::motion_event
```

The system is having trouble getting a good reading due to too much motion.

7.11.2.10 poor_signal

```
bool libct_device_status_t::poor_signal
```

The system failed to calibrate or timed out process signals so measurements were aborted.

7.11.2.11 data_valid

```
bool libct_device_status_t::data_valid
```

There are valid vital signs measurements.

This is used to notify the GUI if data should be displayed or hidden.

7.11.2.12 calibration_offset_failed

```
bool libct_device_status_t::calibration_offset_failed
```

Too much movement.

The calibration offset calculation failed to identify pulses due to movement.

7.11.2.13 weak_signal

```
bool libct_device_status_t::weak_signal
```

Calibration oscillometric curve amplitude is too weak to verify reading.

7.11.2.14 bad_cuff

```
bool libct_device_status_t::bad_cuff
```

The cuff is not holding pressure as expected.

7.11.2.15 ble_adv

```
bool libct_device_status_t::ble_adv
```

The bluetooth module is advertising.

7.11.2.16 recal_soon

```
bool libct_device_status_t::recal_soon
```

An automatic recalibration will be occurring shortly.

7.11.2.17 too_many_fails

```
bool libct_device_status_t::too_many_fails
```

Auto-calibration failed too many consecutive times try manual calibration.

The documentation for this class was generated from the following file:

- caretaker.h

7.12 libct_device_t Class Reference

Handle used to identify a connected device the application is monitoring.

Data Fields

- `int(* get_state)(struct libct_device_t *thiz)`
Return a [device state](#) enumeration representing the current state of the library context that is associated with this device.
- `int(* get_class)(struct libct_device_t *thiz)`
Return the device class that was set in the initialization data passed to [libct_init\(\)](#).
- `const char *(* get_name)(struct libct_device_t *thiz)`
Return the device manufacturer friendly name.
- `const char *(* get_address)(struct libct_device_t *thiz)`
Return the device address.
- `const char *(* get_serial_number)(struct libct_device_t *thiz)`
Return the device serial number.
- `const libct_version_t *(* get_hw_version)(struct libct_device_t *thiz)`
Return the device hardware version.
- `const libct_version_t *(* get_fw_version)(struct libct_device_t *thiz)`
Return the device firmware version.
- `libct_context_t *(* get_context)(struct libct_device_t *thiz)`
Return the library context bound to this device.

7.12.1 Detailed Description

Handle used to identify a connected device the application is monitoring.

The device handle is used to identify and aggregate general information about a connected device, such as the device name, address, serial number, etc., that the application can query. Note each device handle is associated with a library context and you can retrieve it anytime with [libct_get_device\(\)](#) passing the context as argument. As such, you should not hold on to device handles in your application code as they may change when the devices they are associated with become disconnected.

The device handle primary purpose is to identify data notified to your application callbacks.

7.12.2 Field Documentation

7.12.2.1 get_state

```
int ( * libct_device_t::get_state) (struct libct_device_t *thiz)
```

Return a [device state](#) enumeration representing the current state of the library context that is associated with this device.

You would invoke the function as follows, but note a convenience macro with simpler interface than the function is available.

```
// These two calls are similar, but the second is simpler.
// (1) Use the device function.
int state = device->get_state(device);

// (2) Or use macro with simpler interface.
int state = libct_device_get_state(device);
```

Parameters

<i>thiz</i>	The device instance.
-------------	----------------------

See also

[libct_device_get_state\(\)](#)
[libct_device_uninitialized\(\)](#)
[libct_device_intialized\(\)](#)
[libct_device_discovering\(\)](#)
[libct_device_connecting\(\)](#)
[libct_device_connected\(\)](#)
[libct_device_disconnecting\(\)](#)
[libct_device_disconnected\(\)](#)
[libct_device_monitoring\(\)](#)
[libct_device_measuring\(\)](#)

7.12.2.2 get_class

```
int ( * libct_device_t::get_class) (struct libct_device_t *thiz)
```

Return the device class that was set in the initialization data passed to [libct_init\(\)](#).

You would invoke the function as follows, but note a convenience macro with simpler interface than the function is available.

```
// These two calls are similar, but the second is simpler.  
// (1) Use the device function.  
int class = device->get_class(device);  
  
// (2) Or use macro with simpler interface.  
int class = libct_device_get_class(device);
```

Parameters

<i>thiz</i>	The device instance.
-------------	----------------------

See also

[libct_device_get_class\(\)](#)

7.12.2.3 get_name

```
const char* ( * libct_device_t::get_name) (struct libct_device_t *thiz)
```

Return the device manufacturer friendly name.

You would invoke the function as follows, but note a convenience macro with simpler interface than the function is available.

```
// These two calls are similar, but the second is simpler.
// (1) Use the device function.
const char* name = device->get_name(device);

// (2) Or use macro with simpler interface.
const char* name = libct_device_get_name(device);
```

Parameters

<i>thiz</i>	The device instance.
-------------	----------------------

See also

[libct_device_get_name\(\)](#)

7.12.2.4 get_address

```
const char*( * libct_device_t::get_address) (struct libct_device_t *thiz)
```

Return the device address.

You would invoke the function as follows, but note a convenience macro with simpler interface than the function is available.

```
// These two calls are similar, but the second is simpler.
// (1) Use the device function.
const char* address = device->get_address(device);

// (2) Or use macro with simpler interface.
const char* address = libct_device_get_address(device);
```

Parameters

<i>thiz</i>	The device instance.
-------------	----------------------

See also

[libct_device_get_address\(\)](#)

7.12.2.5 get_serial_number

```
const char*( * libct_device_t::get_serial_number) (struct libct_device_t *thiz)
```

Return the device serial number.

You would invoke the function as follows, but note a convenience macro with simpler interface than the function is available.

```
// These two calls are similar, but the second is simpler.  
// (1) Use the device function.  
const char* sn = device->get_serial_number(device);  
  
// (2) Or use macro with simpler interface.  
const char* sn = libct_device_get_serial_number(device);
```

Parameters

<i>thiz</i>	The device instance.
-------------	----------------------

See also

[libct_device_get_serial_number\(\)](#)

7.12.2.6 get_hw_version

```
const libct_version_t*( * libct_device_t::get_hw_version) (struct libct_device_t *thiz)
```

Return the device hardware version.

You would invoke the function as follows, but note a convenience macro with simpler interface than the function is available.

```
// These two calls are similar, but the second is simpler.
// (1) Use the device function.
const libct_version_t* version = device->get_hw_version(device);

// (2) Or use macro with simpler interface.
const libct_version_t* version = libct_device_get_hw_version(
    device);
```

Parameters

<i>thiz</i>	The device instance.
-------------	----------------------

See also

[libct_device_get_hw_version\(\)](#)

7.12.2.7 get_fw_version

```
const libct_version_t*( * libct_device_t::get_fw_version) (struct libct_device_t *thiz)
```

Return the device firmware version.

You would invoke the function as follows, but note a convenience macro with simpler interface than the function is available.

```
// These two calls are similar, but the second is simpler.
// (1) Use the device function.
const libct_version_t* version = device->get_fw_version(device);

// (2) Or use macro with simpler interface.
const libct_version_t* version = libct_device_get_fw_version(
    device);
```

Parameters

<i>thiz</i>	The device instance.
-------------	----------------------

See also

[libct_device_get_fw_version\(\)](#)

7.12.2.8 get_context

```
libct_context_t* ( * libct_device_t::get_context ) (struct libct_device_t *thiz)
```

Return the library context bound to this device.

You would invoke the function as follows, but note a convenience macro with simpler interface than the function is available.

```
// These two calls are similar, but the second is simpler.
// (1) Use the device function.
libct_context_t* context = device->get_context(device);

// (2) Or use macro with simpler interface.
libct_context_t* context = libct_device_get_context(device);
```

Parameters

<i>thiz</i>	The device instance.
-------------	----------------------

The documentation for this class was generated from the following file:

- caretaker.h

7.13 libct_init_data_t Struct Reference

Structure defining initialization data passed to [libct_init\(\)](#).

Data Fields

- int [device_class](#)
The device class.

7.13.1 Detailed Description

Structure defining initialization data passed to [libct_init\(\)](#).

7.13.2 Field Documentation

7.13.2.1 device_class

```
int libct_init_data_t::device_class
```

The [device class](#).

The documentation for this struct was generated from the following file:

- caretaker.h

7.14 libct_monitor_flags_t Class Reference

Data monitor flags passed to [libct_start_monitoring\(\)](#)

7.14.1 Detailed Description

Data monitor flags passed to [libct_start_monitoring\(\)](#)

The documentation for this class was generated from the following file:

- caretaker.h

7.15 libct_param_pulse_t Class Reference

Parametrized pulse data within the [libct_stream_data_t](#) packet.

Data Fields

- bool [valid](#)
The other fields are valid when this field is non-zero (true) and invalid otherwise.
- short [protocol_header](#)
Date transfer control byte used internally to assemble the data.
- short [t0](#)
Pulse onset time (index).
- short [t1](#)
First pulse peak time (index).
- short [t2](#)
Second pulse peak time (index).
- short [t3](#)
Third pulse peak time (index).
- int [p0](#)

- Integrated pulse onset value.*
- int [p1](#)
First integrated pulse peak value.
- int [p2](#)
Second integrated pulse peak value.
- int [p3](#)
Third integrated pulse peak value.
- short [ibi](#)
Inter-beat interval (1/HR) in samples @ 500Hz.
- short [as](#)
Arterial stiffness.
- short [sqe](#)
Signal quality estimate.
- short [pressure](#)
The most recent raw ADC cuff pressure.
- unsigned int [time](#)
Relative system time of occurrence.
- int [waveform_len](#)
The number of signed int8 snapshot data points.
- char [waveform](#) [0]
The pulse snapshot waveform data.

7.15.1 Detailed Description

Parametrized pulse data within the [libct_stream_data_t](#) packet.

The parametrized pulse data is an aggregate of the pulse parameters and pulse snapshot waveform data.

7.15.2 Field Documentation

7.15.2.1 valid

```
bool libct_param_pulse_t::valid
```

The other fields are valid when this field is non-zero (true) and invalid otherwise.

7.15.2.2 protocol_header

```
short libct_param_pulse_t::protocol_header
```

Date transfer control byte used internally to assemble the data.

7.15.2.3 t0

```
short libct_param_pulse_t::t0
```

Pulse onset time (index).

7.15.2.4 t1

```
short libct_param_pulse_t::t1
```

First pulse peak time (index).

7.15.2.5 t2

```
short libct_param_pulse_t::t2
```

Second pulse peak time (index).

7.15.2.6 t3

```
short libct_param_pulse_t::t3
```

Third pulse peak time (index).

7.15.2.7 p0

```
int libct_param_pulse_t::p0
```

Integrated pulse onset value.

7.15.2.8 p1

```
int libct_param_pulse_t::p1
```

First integrated pulse peak value.

7.15.2.9 p2

```
int libct_param_pulse_t::p2
```

Second integrated pulse peak value.

7.15.2.10 p3

```
int libct_param_pulse_t::p3
```

Third integrated pulse peak value.

7.15.2.11 ibi

```
short libct_param_pulse_t::ibi
```

Inter-beat interval (1/HR) in samples @ 500Hz.

7.15.2.12 as

```
short libct_param_pulse_t::as
```

Arterial stiffness.

7.15.2.13 sqe

```
short libct_param_pulse_t::sqe
```

Signal quality estimate.

7.15.2.14 pressure

```
short libct_param_pulse_t::pressure
```

The most recent raw ADC cuff pressure.

7.15.2.15 time

```
unsigned int libct_param_pulse_t::time
```

Relative system time of occurrence.

7.15.2.16 waveform_len

```
int libct_param_pulse_t::waveform_len
```

The number of signed int8 snapshot data points.

The documentation for this class was generated from the following file:

- caretaker.h

7.16 libct_posture_t Class Reference

Patient postures.

7.16.1 Detailed Description

Patient postures.

The documentation for this class was generated from the following file:

- caretaker.h

7.17 libct_pulse_ox_t Class Reference

Pulse oximetry data point within the [libct_stream_data_t](#) packet.

Data Fields

- bool [valid](#)
The other fields are valid when this field is non-zero (true) and invalid otherwise.
- int [sao2](#)
Blood oxygen level (percentage).
- int [pulse_rate](#)
Pulse rate in beats per minute (30-200BPM).
- unsigned int [timestamp](#)
Time stamp from the device associated with the data.

7.17.1 Detailed Description

Pulse oximetry data point within the [libct_stream_data_t](#) packet.

Note

Reserved for future use.

7.17.2 Field Documentation

7.17.2.1 sao2

```
int libct_pulse_ox_t::sao2
```

Blood oxygen level (percentage).

7.17.2.2 pulse_rate

```
int libct_pulse_ox_t::pulse_rate
```

Pulse rate in beats per minute (30-200BPM).

7.17.2.3 timestamp

```
unsigned int libct_pulse_ox_t::timestamp
```

Time stamp from the device associated with the data.

The documentation for this class was generated from the following file:

- caretaker.h

7.18 libct_pulse_t Class Reference

Raw or integrated pulse data point within the [libct_stream_data_t](#) packet.

Data Fields

- bool [valid](#)
The other fields are valid when this field is non-zero (true) and invalid otherwise.
- short [value](#)
Pulse value.
- unsigned int [timestamp](#)
Counter value associated with the pulse value.

7.18.1 Detailed Description

Raw or integrated pulse data point within the [libct_stream_data_t](#) packet.

7.18.2 Field Documentation

7.18.2.1 timestamp

```
unsigned int libct_pulse_t::timestamp
```

Counter value associated with the pulse value.

The documentation for this class was generated from the following file:

- caretaker.h

7.19 libct_pulse_waveform_t Struct Reference

Pulse waveform data returned from the device as a result of a previous read request.

Data Fields

- ```
struct {
 libct_pulse_t * datapoints
 unsigned int count
} int_pulse
```

  
*Array of integrated waveform pulse data points.*
- ```
struct {  
    libct\_param\_pulse\_t * datapoints  
    unsigned int count  
} param\_pulse
```


Array of parameterize pulse waveform data points.
- long [receive_time](#)
Value of local clock (in milliseconds) when this stream packet was received.

7.19.1 Detailed Description

Pulse waveform data returned from the device as a result of a previous read request.

Note this data is returned only after an explicit read of raw pulse waveform data.

7.19.2 Field Documentation

7.19.2.1 receive_time

```
long libct_pulse_waveform_t::receive_time
```

Value of local clock (in milliseconds) when this stream packet was received.

This time stamp is used to measure processing latency and for history logging. It differs from time stamp found in each data point generated by the remote device.

The documentation for this struct was generated from the following file:

- caretaker.h

7.20 libct_status_t Class Reference

Function return status codes.

7.20.1 Detailed Description

Function return status codes.

The documentation for this class was generated from the following file:

- caretaker.h

7.21 libct_stream_data_t Class Reference

This structure is used to hand-off data received from the remote device to the application.

Data Fields

- [libct_device_t](#) * [device](#)
Reference to the device that generated this data.
- [libct_device_status_t](#) [device_status](#)
Device status information.
- [libct_battery_info_t](#) [battery_info](#)
Battery information.
- struct {
 [libct_vitals_t](#) * [datapoints](#)
 Array of vital sign datapoints.
 unsigned int [count](#)
 The count of datapoints.
 } [vitals](#)

Array of vital sign data points.
- struct {
 [libct_cuff_pressure_t](#) * [datapoints](#)
 Array of cuff pressure data points.
 unsigned int [count](#)
 The count of data points.
 } [cuff_pressure](#)

Array of cuff pressure data points.
- struct {
 [libct_temperature_t](#) * [datapoints](#)
 Array of temperature data points.
 unsigned int [count](#)
 The count of data points.
 } [temperature](#)

Array of temperature data points.
- struct {
 [libct_pulse_ox_t](#) * [datapoints](#)
 Array of spo2 data points.
 unsigned int [count](#)
 The count of data points.
 } [pulse_ox](#)

Array of pulse oximetry data points.
- struct {
 [libct_vitals2_t](#) * [datapoints](#)
 Array of secondary vital sign data points.
 unsigned int [count](#)
 The count of data points.
 } [vitals2](#)

Array of secondary vital sign data points.
- struct {
 [libct_pulse_t](#) * [datapoints](#)
 Array of raw pulse (pulse rate) data points.
 unsigned int [count](#)
 The count of data points.
 } [raw_pulse](#)

Array of raw pulse (pulse rate) waveform data points.

- struct {
 - [libct_pulse_t](#) * [datapoints](#)
Array of integrated pulse (pulse pressure) data points.
 - unsigned int [count](#)
The count of data points.
- } [int_pulse](#)
- Array of integrated pulse (pulse pressure) waveform data points.
- struct {
 - [libct_param_pulse_t](#) * [datapoints](#)
Array of parameterize pulse (pulse snapshot) data points.
 - unsigned int [count](#)
The count of data points.
- } [param_pulse](#)
- Array of parameterized pulse data.
- struct {
 - [libct_cal_curve_t](#) * [datapoints](#)
Array of calibration curve data points.
 - unsigned int [count](#)
The count of data points.
- } [cal_curve](#)
- Array of calibration curve data points.
- long [receive_time](#)
Value of local clock (in milliseconds) when this stream packet was received.

7.21.1 Detailed Description

This structure is used to hand-off data received from the remote device to the application.

Data from the device is sent automatically after calling [libct_start_monitoring\(\)](#) successfully, and delivered to your application via the [on_data_received\(\)](#) callback function. This data structure is a container of arrays grouping one or more records of the same data type at different time instances. The various array data types are not produced coherently at the device so not all fields will be populated in stream data packets delivered to the application. If no data is available for a given array, the array data points field will be set to null and the count set to zero to signal no data.

The stream data packets notified to the application depends on the [monitor flags](#) passed to [libct_start_monitoring\(\)](#) and whether or not [libct_start_measuring\(\)](#) was called to start taking vital sign measurements. So you can control the data reported to the application by specifying only the monitoring flags corresponding to the data you care about.

With the exception of the [device_status](#) and [battery_info](#) data members that are not array fields, the following convenience macros are available to access array entries within the stream data packet. More details about usage is provided in the description for each stream data member where the macros apply.

- [libct_dp_count\(\)](#)
- [libct_get_dp\(\)](#)
- [libct_get_first_dp\(\)](#)
- [libct_get_last_dp\(\)](#)
- [for_each_dp\(\)](#)

7.21.2 Field Documentation

7.21.2.1 device

```
libct_device_t* libct_stream_data_t::device
```

Reference to the device that generated this data.

7.21.2.2 device_status

```
libct_device_status_t libct_stream_data_t::device_status
```

Device status information.

7.21.2.3 battery_info

```
libct_battery_info_t libct_stream_data_t::battery_info
```

Battery information.

7.21.2.4 datapoints [1/8]

```
libct_vitals_t* libct_stream_data_t::datapoints
```

Array of vital sign datapoints.

7.21.2.5 count

```
unsigned int libct_stream_data_t::count
```

The count of datapoints.

The count of data points.

7.21.2.6 vitals

```
struct { ... } libct_stream_data_t::vitals
```

Array of vital sign data points.

For convenience, you can use the macros [libct_get_last_dp\(\)](#), [libct_get_first_dp\(\)](#), and [libct_get_dp\(\)](#) to extract a single vital sign data point from the stream packet like so.

```
libct_vitals_t* dp = libct_get_last_dp(data, vitals);
if ( dp && dp->valid ) {
    // use most recent vitals data point
}
```

Alternatively, you could iterate over all vital sign data points with the [for_each_dp\(\)](#) macro like so.

```
libct_vitals_t* dp;
unsigned int idx;
for_each_dp(data, idx, dp, vitals) {
    if ( dp && dp->valid ) {
        // use vitals data point
    }
}
```

7.21.2.7 datapoints [2/8]

```
libct_cuff_pressure_t* libct_stream_data_t::datapoints
```

Array of cuff pressure data points.

7.21.2.8 cuff_pressure

```
struct { ... } libct_stream_data_t::cuff_pressure
```

Array of cuff pressure data points.

For convenience, you can use the macros [libct_get_last_dp\(\)](#), [libct_get_first_dp\(\)](#), and [libct_get_dp\(\)](#) to extract a single cuff pressure data point from the stream packet like so.

```
libct_cuff_pressure_t* dp = libct_get_last_dp(data,
    cuff_pressure);
if ( dp && dp->valid ) {
    // use most recent cuff pressure data point
}
```

Alternatively, you could iterate over all cuff pressure data points with the [for_each_dp\(\)](#) macro like so.

```
libct_cuff_pressure_t* dp;
unsigned int idx;
for_each_dp(data, idx, dp, cuff_pressure) {
    if ( dp && dp->valid ) {
        // use cuff pressure data point
    }
}
```

7.21.2.9 datapoints [3/8]

```
libct_temperature_t* libct_stream_data_t::datapoints
```

Array of temperature data points.

7.21.2.10 temperature

```
struct { ... } libct_stream_data_t::temperature
```

Array of temperature data points.

For convenience, you can use the macros [libct_get_last_dp\(\)](#), [libct_get_first_dp\(\)](#), and [libct_get_dp\(\)](#) to extract a single temperature data point from the stream packet like so.

```
libct_temperature_t* dp = libct_get_last_dp(data,
    temperature);
if ( dp && dp->valid ) {
    // use most recent temperature data point
}
```

Alternatively, you could iterate over all temperature data points with the [for_each_dp\(\)](#) macro like so.

```
libct_temperature_t* dp;
unsigned int idx;
for_each_dp(data, idx, dp, temperature) {
    if ( dp && dp->valid ) {
        // use temperature data point
    }
}
```

7.21.2.11 datapoints [4/8]

```
libct_pulse_ox_t* libct_stream_data_t::datapoints
```

Array of spo2 data points.

7.21.2.12 pulse_ox

```
struct { ... } libct_stream_data_t::pulse_ox
```

Array of pulse oximetry data points.

For convenience, you can use the macros [libct_get_last_dp\(\)](#), [libct_get_first_dp\(\)](#), and [libct_get_dp\(\)](#) to extract a single spo2 data point from the stream packet like so.

```
libct_pulse_ox_t* dp = libct_get_last_dp(data,
    pulse_ox);
if ( dp && dp->valid ) {
    // use most recent spo2 data point
}
```

Alternatively, you could iterate over all spo2 data points with the [for_each_dp\(\)](#) macro like so.

```
libct_pulse_ox_t* dp;
unsigned int idx;
for_each_dp(data, idx, dp, pulse_ox) {
    if ( dp && dp->valid ) {
        // use spo2 data point
    }
}
```

7.21.2.13 datapoints [5/8]

```
libct_vitals2_t* libct_stream_data_t::datapoints
```

Array of secondary vital sign data points.

7.21.2.14 vitals2

```
struct { ... } libct_stream_data_t::vitals2
```

Array of secondary vital sign data points.

Note

The secondary vitals are for internal use or research only.

For convenience, you can use the macros [libct_get_last_dp\(\)](#), [libct_get_first_dp\(\)](#), and [libct_get_dp\(\)](#) to extract a single secondary vital sign data point from the stream packet.

```
libct_vitals2_t* dp = libct_get_last_dp(data,
    vitals2);
if ( dp && dp->valid ) {
    // use most recent secondary vitals data point
}
```

Alternatively, you could iterate over all secondary vital sign data points with the [for_each_dp\(\)](#) macro like so.

```
libct_vitals2_t* dp;
unsigned int idx;
for_each_dp(data, idx, dp, vitals) {
    if ( dp && dp->valid ) {
        // use secondary vitals data point
    }
}
```

7.21.2.15 datapoints [6/8]

```
libct_pulse_t* libct_stream_data_t::datapoints
```

Array of raw pulse (pulse rate) data points.

Array of integrated pulse (pulse pressure) data points.

7.21.2.16 raw_pulse

```
struct { ... } libct_stream_data_t::raw_pulse
```

Array of raw pulse (pulse rate) waveform data points.

For convenience, you could iterate over all raw pulse data points with the [for_each_dp\(\)](#) macro like so.

```
libct_pulse_t* dp;
unsigned int idx;
for_each_dp(data, idx, dp, raw_pulse) {
    if ( dp && dp->valid ) {
        // use pulse data point
    }
}
```

7.21.2.17 int_pulse

```
struct { ... } libct_stream_data_t::int_pulse
```

Array of integrated pulse (pulse pressure) waveform data points.

For convenience, you could iterate over all integrated pulse data points with the [for_each_dp\(\)](#) macro like so.

```
libct_pulse_t* dp;
unsigned int idx;
for_each_dp(data, idx, dp, int_pulse) {
    if ( dp && dp->valid ) {
        // use pulse data point
    }
}
```

7.21.2.18 datapoints [7/8]

```
libct_param_pulse_t* libct_stream_data_t::datapoints
```

Array of parameterize pulse (pulse snapshot) data points.

7.21.2.19 param_pulse

```
struct { ... } libct_stream_data_t::param_pulse
```

Array of parameterized pulse data.

The data is an aggregate of the pulse parameters and pulse snapshot waveform data.

For convenience, you can use the macros [libct_get_last_dp\(\)](#), [libct_get_first_dp\(\)](#), and [libct_get_dp\(\)](#) to extract a single pulse snapshot from the stream packet like so.

```
libct_param_pulse_t* dp = libct_get_last_dp(data,
    param_pulse);
if ( dp && dp->valid ) {
    // use most recent pulse snapshot
}
```

Alternatively, you could iterate over all pulse snapshot data points with the [for_each_dp\(\)](#) macro like so.

```
libct_param_pulse_t* dp;
unsigned int idx;
for_each_dp(data, idx, dp, pulse_param) {
    if ( dp && dp->valid ) {
        // use pulse snapshot
    }
}
```

7.21.2.20 datapoints [8/8]

```
libct_cal_curve_t* libct_stream_data_t::datapoints
```

Array of calibration curve data points.

7.21.2.21 cal_curve

```
struct { ... } libct_stream_data_t::cal_curve
```

Array of calibration curve data points.

Note

The cal curve data is for internal use or research only

For convenience, you can use the macros [libct_get_last_dp\(\)](#), [libct_get_first_dp\(\)](#), and [libct_get_dp\(\)](#) to extract a single calibration curve data point from the stream packet like so.

```
libct_cal_curve_t* dp = libct_get_last_dp(data,
    cal_curve);
if ( dp && dp->valid ) {
    // use most recent calibration curve data point
}
```

Alternatively, you could iterate over all calibration curve data points with the [for_each_dp\(\)](#) macro like so.

```
libct_cal_curve_t* dp;
unsigned int idx;
for_each_dp(data, idx, dp, cal_curve) {
    if ( dp && dp->valid ) {
        // use calibration curve data point
    }
}
```

7.21.2.22 receive_time

```
long libct_stream_data_t::receive_time
```

Value of local clock (in milliseconds) when this stream packet was received.

This time stamp is used to measure processing latency and for history logging. It differs from time stamp found in each data point generated by the remote device.

The documentation for this class was generated from the following file:

- caretaker.h

7.22 libct_temperature_t Class Reference

Temperature data point within the [libct_stream_data_t](#) packet.

Data Fields

- bool [valid](#)
The other fields are valid when this field is non-zero (true) and invalid otherwise.
- int [value](#)
The temperature value.
- unsigned int [timestamp](#)
Time stamp from the device associated with the data.

7.22.1 Detailed Description

Temperature data point within the [libct_stream_data_t](#) packet.

Note

Reserved for future use.

7.22.2 Field Documentation

7.22.2.1 value

```
int libct_temperature_t::value
```

The temperature value.

7.22.2.2 timestamp

```
unsigned int libct_temperature_t::timestamp
```

Time stamp from the device associated with the data.

The documentation for this class was generated from the following file:

- caretaker.h

7.23 libct_version_t Struct Reference

CareTaker version information.

Data Fields

- int [major](#)
major version number
- int [minor](#)
minor version number
- int [revision](#)
revision number
- int [build](#)
build number

7.23.1 Detailed Description

CareTaker version information.

The documentation for this struct was generated from the following file:

- caretaker.h

7.24 libct_vitals2_t Class Reference

Secondary Vitals data point within the [libct_stream_data_t](#) packet.

Data Fields

- bool [valid](#)
The other fields are valid when this field is non-zero (true) and invalid otherwise.
- unsigned short [blood_volume](#)
Blood volume in mS.
- unsigned char [cardiac_output](#)
Cardiac output in L/min.
- unsigned short [ibi](#)
Inter-beat Interval in mS.
- unsigned short [lvet](#)
Left ventricular ejection time.
- float [p2p1](#)
P ratio.
- float [hrComp](#)
hrComp
- float [pr](#)
pr
- int [reserved](#) [7]
Reserved for future use.
- unsigned int [timestamp](#)
Time stamp from the device associated with the data.

7.24.1 Detailed Description

Secondary Vitals data point within the [libct_stream_data_t](#) packet.

Note

The secondary vitals are for internal use or research only.

7.24.2 Field Documentation

7.24.2.1 blood_volume

```
unsigned short libct_vitals2_t::blood_volume
```

Blood volume in mS.

7.24.2.2 timestamp

```
unsigned int libct_vitals2_t::timestamp
```

Time stamp from the device associated with the data.

The documentation for this class was generated from the following file:

- caretaker.h

7.25 libct_vitals_t Class Reference

Vitals data point within the [libct_stream_data_t](#) packet.

Data Fields

- bool [valid](#)
The other fields are valid when this field is non-zero (true) and invalid otherwise.
- bool [bp_status](#)
An indicator if a valid blood pressure was found or if the algorithm failed.
- bool [map_status](#)
An indicator of if a valid MAP measurement has been integrated.
- bool [hr_status](#)
An indicator if a valid HR has been determined.
- bool [respiration_status](#)
An indicator if a valid respiration reading was found.
- bool [integration_error](#)
General catchall for integration errors.
- bool [differentiation_error](#)
A discontinuity was detected in the differentiation.
- bool [p12_finder_error](#)
Unable to locate P1 P2 within the pulse.
- bool [p3_finder_eError](#)
Unable to locate P3 within the pulse.
- bool [min_index_out_of_range](#)
The onset of the pulse was not found in the allowable window, so the values are being discarded.
- bool [max_index_out_of_range](#)
The index of the minimum point in the integral was out of range.
- bool [slope_out_of_range](#)
The slope correction of the integrated pulse was out of range.
- short [systolic](#)
Systolic measurement.
- short [diastolic](#)
Diastolic measurement.
- short [map](#)
Mean arterial pressure value.
- short [heart_rate](#)
Heart rate measurement.
- short [respiration](#)
Respiration measurement.
- short [as](#)
AS factor.
- short [sqe](#)
Signal quality estimate (sqe).
- unsigned int [timestamp](#)
Time stamp from the device associated with the data.

7.25.1 Detailed Description

Vitals data point within the [libct_stream_data_t](#) packet.

7.25.2 Field Documentation

7.25.2.1 valid

```
bool libct_vitals_t::valid
```

The other fields are valid when this field is non-zero (true) and invalid otherwise.

7.25.2.2 bp_status

```
bool libct_vitals_t::bp_status
```

An indicator if a valid blood pressure was found or if the algorithm failed.

True indicates pulse information is valid.

7.25.2.3 systolic

```
short libct_vitals_t::systolic
```

Systolic measurement.

7.25.2.4 diastolic

```
short libct_vitals_t::diastolic
```

Diastolic measurement.

7.25.2.5 map

```
short libct_vitals_t::map
```

Mean arterial pressure value.

7.25.2.6 heart_rate

```
short libct_vitals_t::heart_rate
```

Heart rate measurement.

7.25.2.7 respiration

```
short libct_vitals_t::respiration
```

Respiration measurement.

7.25.2.8 as

```
short libct_vitals_t::as
```

AS factor.

7.25.2.9 sqe

```
short libct_vitals_t::sqe
```

Signal quality estimate (sqe).

Values are in the range [0, 1000], so the sqe can be expressed relatively as a percentage by dividing by 10, .i.e. sqe/10 %.

7.25.2.10 timestamp

```
unsigned int libct_vitals_t::timestamp
```

Time stamp from the device associated with the data.

The documentation for this class was generated from the following file:

- caretaker.h